
BAIRD
Max 30™

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Instructions for Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Description

The **Bard® Max 30™** Inflation Device is a one-piece, plastic, 20cc disposable inflation device with a lock lever design that controls the piston, a manometer, and a high-pressure connecting tube with a male rotating adapter. Also enclosed is a 3-way stopcock for use during preparation of the device. The manometer measures pressures ranging from vacuum to 30 atm; the gauge is marked in 1 atm increments. The accuracy of the manometer has been determined to be within 1 atm over the range.

How Supplied

The **Bard Max 30** inflation device is supplied sterile unless the package has been opened or damaged. For single use only. Do not reuse. Do not resterilize.

Indications for Use

The **Bard Max 30** inflation device is recommended for use while performing peripheral balloon angioplasty to inflate the balloon, monitor the pressure within the balloon and deflate the balloon.

Contraindications

None known.

Warnings

- Always follow the manufacturer's directions accompanying the balloon dilatation catheter for instruction for use, maximum balloon inflation pressure, precautions, and warnings for that device.**
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.**

Precautions

Before use, inspect the device to verify that no damage has occurred during shipping and handling.

Directions for Use**Handling and Storage**

Store in a cool, dry, dark place.

Equipment Required

- Contrast medium
- Sterile saline solution

Preparation

Make all aspiration and injection maneuvers with the lock lever pushed left, i.e., unlocked.

Unlock the piston by pushing the lock lever left. In this position, you can freely pull the piston back for aspiration, or push it forward for injection. To lock the piston in position, slide the lever right to the straight up position.

- Prepare a solution of contrast medium and normal saline in a small sterile bowl or in the well provided with the package. Check balloon catheter instructions for specific contrast mixture recommendations.
- Orient the tubing downward into the contrast medium.
- Push the release lever left and aspirate enough solution to fill the syringe. Attach stopcock.
- Hold the device upright to purge the air from the syringe and connecting tube. Tap the syringe lightly, if necessary, to remove all the air bubbles and fill the connecting tube completely.

- Inspect the syringe, tubing and stopcock to ensure that the device has been completely purged of air bubbles.
- Adjust the syringe volume to the desired amount. If more contrast solution is needed, submerge it into the basin of solution and aspirate. Close stopcock.

Attaching the inflation device to the peripheral balloon dilatation catheter

- Prepare and test the peripheral balloon dilatation catheter according to the manufacturer's directions for use.
- If a separate syringe was used to prepare the balloon catheter, remove it. Create a fluid-fluid connection between the catheter connector and the connecting tube of the inflation device by opening stopcock and placing a drop of contrast solution from the syringe into each hub.
- Hand tighten the hubs securely.

Balloon inflation and deflation

- Release the lock lever and allow the piston to move forward into neutral position (0 atm).
- To inflate the balloon, engage the lock lever, turn the palm grip clockwise slowly until the desired inflation pressure is reached. The lock lever maintains the increasing pressure.
- To deflate the balloon, push the lock lever left, releasing the piston, and pull back. Slide the lock lever back to lock, if desired.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Copyright © 2007 C. R. Bard, Inc. All rights reserved. Printed in the USA.
Copyright © 2007 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2007 C. R. Bard, Inc. Alle Rechte vorbehalten. Druck: USA.
Copyright © 2007 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2007 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright © 2007, C.R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.