Instructions for Use

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PERIPHERAL
VASCULAR

INSTRUCTIONS FOR USE

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

Venaflo® ePTFE Vascular Graft

Venaflo Graft with Carbon

Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

Device Description

The Venaflo ePTFE vascular graft is an IMPRA® expanded polytetrafluoroethylene (ePTFE) vascular graft with a modified venous end. The Venaflo graft with carbon has carbon impregnated into the inner portions of the graft wall.

All IMPRA ePTFE vascular grafts are supplied sterile and non-pyrogenic unless the package is opened or damaged. IMPRA ePTFE vascular grafts are sterilized by ethylene oxide. Each graft is intended for single patient use only. **DO NOT RESTERILIZE.**

Indications for Use

Venaflo ePTFE vascular grafts in Straight, Stepped, CenterFlex™, and Stepped CenterFlex graft configurations are indicated for use as a subcutaneous arteriovenous conduit for blood access only. Stepped configurations may help minimize the risk of steal syndrome and high cardiac output. CenterFlex graft configurations have an external spiral bead support on the center of the graft and can be used where resistance to compression or kinking is desired.

Contraindications

None known.

Note: Numbers in the "Warnings" and "Precautions" sections are for ease of reference only.

Warnings

1. Do NOT cannulate the externally supported portion of any Venaflo ePTFE vascular graft due to the risk of beading embolization and/or pseudoaneurysm formation.
2. DO NOT REMOVE THE EXTERNAL SPIRAL BEAD SUPPORT FROM ANY CENTERFLEX CONFIGURED GRAFT. Attempts to remove the bead may damage the graft. If damage occurs, discard the graft.
3. ePTFE vascular grafts are not elastic in the longitudinal direction. Failure to cut the grafts to an appropriate length may result in anastomotic disruption or excessive bleeding.
4. Do NOT expose IMPRA ePTFE products to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.¹
5. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

Precautions

1. Venaflo ePTFE vascular grafts are not radially elastic. When embolectomy or balloon angioplasty catheters are used within the lumen of a Venaflo graft, match the inflated balloon size to the inner diameter of the graft. Ensure the graftotomy is large enough to extract a fully expanded balloon. Over-inflation of the balloon or use of an inappropriately sized balloon catheter may dilate or damage the graft.
2. If clamping is necessary, use only atraumatic or appropriate vascular smooth jawed clamps to avoid damage to the graft wall. Avoid repeated or excessive clamping at the same location on the graft. Do not clamp the cuffed portion of the graft.
3. Preclotting of this graft is unnecessary. Exposure to any solution (alcohol, oil, aqueous solutions, etc.) prior to implantation may result in seroma formation when the graft is implanted.²
4. Do not forcibly inject any solution through the lumen of the graft, or fill the graft with fluid prior to pulling it through the tunnel as this may result in seroma formation. For more information, ask for IMPRA Technical Report TR-100.
5. When suturing, avoid excessive tension on the suture line, gaps between the graft and host vessel, and inappropriate suture spacing and bites. Failure to avoid the above items may lead to anastomotic bleeding. For further instructions, refer to the "Suturing" section.
6. Do not cannulate the cuffed portion of the graft. For further instructions, refer to the "Blood Access" section.

Adverse Reactions:

Potential complications which may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: disruption or tearing of the suture line, graft, and/or host vessel; suture hole bleeding; thrombosis; seroma formation; swelling of the implanted limb; formation of hematomas or pseudoaneurysms; infection; steal syndrome; and/or skin erosion.

DIRECTIONS FOR USE

Equipment Required:

Tunneler, double-armed suture, atraumatic clamp, scalpel and/or scissors.

Opening the Package:

While holding the outer tray in one hand, peel back the lid and remove the inner tray. Peel back the inner tray lid slowly and carefully remove the graft using sterile atraumatic instruments or gloves. Protect the graft against damage from sharp or heavy instruments.

General Operative Techniques


Venous Anastomotic Preparation:

Use the following guidelines when sizing the cuffed end of all Venaflo ePTFE vascular grafts for the venous anastomosis.

Size of Vein (external diameter)	Cuff Sizing (see Figure 1)
> 8 mm	Do not trim
6-8 mm	Trim 5 mm off toe
3-5 mm	Trim 10 mm off toe

Trimming Recommendations

Flatten the cuff between the thumb and finger to obtain a cuff profile. Using curved scissors, trim the cuff with the curvature suggested by the dotted lines in Figure 1. Start at the heel when trimming, being careful not to trim any material from the heel, as this is needed for suturing.

For more information, ask for IMPRA Technical Report TR-202.

Operative Techniques (Use of Kelly-Wick Tunneler or AV Access Sheath Tunneler)
(Refer also to the Kelly-Wick or AV Access Sheath Tunneler Directions for Use)

Tunneling Recommendations

Venaflo ePTFE Vascular Grafts in Straight Anatomical Configurations:

Due to the larger area of the Venaflo graft's cuff, it may be useful to perform the venous anastomosis before tunneling.

After selecting the appropriate tunneler, introduce the tunneler through the arterial site. After creating the tunnel, feed the graft from the venous end to the arterial end OR attach the arterial end of the graft to the tunneler rod and pull the graft through the tunnel (see Figure 2). Ensure that the cuffed end does not pass through the tunnel. Only the arterial end of the Venaflo ePTFE vascular graft should be sized for length.

Figure 2

Kelly-Wick Bullet-tip Selection: Use the following guidelines to select a bullet tip.

1. For Straight (non-flexed) grafts, select a bullet tip that matches the internal diameter of the graft to create the tunnel. Use the same bullet tip to pull the graft through the tunnel, attaching it to the arterial (non-cuffed) graft end.
2. For Stepped grafts, use a bullet tip that matches the largest stated internal diameter of the graft to create the tunnel, then switch to a bullet tip that matches the arterial (non-cuffed) end of the graft to pull it through the tunnel.
3. For CenterFlex and Stepped CenterFlex grafts, use a bullet tip one or two millimeters greater than the largest stated internal diameter of the graft to create the tunnel. After creating the tunnel, replace the bullet tip with a tip of the same size as the internal diameter of the arterial (non-cuffed) end of the graft. Attach the arterial end of the graft to the tunneler with this tip and pull the graft through the tunnel.

Venaflor® ePTFE Vascular Grafts in Looped Anatomical Configurations:

Due to the larger area of the **Venaflor** graft's cuff, it may be useful to perform the venous anastomosis before tunneling.

Perform an incision to expose the artery and vein to be anastomosed and a counter-incision at the site where the graft will curve. For **CenterFlex™** grafts, this is where the flexed portion of the graft should be placed. Position the graft overlying the forearm to optimize graft placement prior to tunneling. Introduce the tunneler at the counter-incision site. After creating the first tunnel, replace the bullet tip (if necessary), attach the arterial end of the graft to the tunneler, and pull the graft through the tunnel leaving the cuffed end exposed near the vein to be anastomosed.

For the second tunnel, repeat the tunneling procedure, initiating the tunnel at the anastomotic site and advancing it down to the counter-incision (see **Figure 3**). Once again, ensure that the venous (cuffed) graft end does not pass into the tunnel.

Figure 3

For more information, ask for **IMPRA®** Technical Report TR-201 or TR-203.

Suturing

Best results are achieved by using a tapered, noncutting needle with a nonabsorbable monofilament suture approximately the same size as the needle. Take 1 to 2 mm suture bites in the graft.

Follow the curve of the needle and pull the suture at a 90° angle to minimize suture hole elongation and bleeding. Proper sizing of the graft prior to implant will minimize suture hole elongation caused by excessive tension. Application of a topical hemostatic agent may aid in controlling any bleeding that may occur.

See **PRECAUTION #5**.

Suggested Suturing Technique

After sizing the cuff to the vein diameter (refer to the "Venous Anastomotic Preparation" section), perform a venotomy five to eight millimeters shorter than the length of the cuff.

Using double armed sutures, tack the heel of the cuff, aligning the toe in the direction of the venous blood flow.

Suture continuously from the heel to the midpoint of the cuff on one side of the venotomy. Using the other arm of the heel suture, suture to the midpoint of the graft on the other side of the venotomy.

Extend the venotomy to the exact cuff length (see **Figure 4**).

Figure 4

Complete the "cuffed" anastomosis by suturing with one of the suture arms, around the toe, to the other side of the cuff. Always tie off on the side of the cuff (see **Figure 5**).

Figure 5

Perform the arterial anastomosis using standard surgical techniques. In order to retain a cuffed configuration, only the arterial end of the **Venaflor** ePTFE vascular grafts should be sized for length. See **WARNING #3**.

For more information, ask for **IMPRA** Technical Report TR-202.

Thrombectomy

In the event of a postoperative occlusion, **IMPRA** grafts can be declotted (see **PRECAUTION #1**) using either a longitudinal or transverse incision. During the early postoperative period, the natural progression of wound healing renders the graft translucent in appearance. When the material is in this state, it is recommended that a longitudinal incision with stay sutures be used for thrombectomy procedures. If a transverse incision is performed, the use of a horizontal mattress suture technique may be required to aid in closure.

Longitudinal Incision:

Place stay sutures before introducing the embolectomy catheter. (If the graft is a **CenterFlex** product, cut through the spiral beading and base tube. After closure, the spiral beading will realign itself.) Place a longitudinal incision in the graft that is long enough to accommodate the extraction of a fully dilated thrombectomy catheter balloon.

Transverse Incision:

No stay sutures are necessary and a horizontal mattress suture technique will aid in closure.

Angiography

Should angiography be performed at the time of the procedure, the artery proximal to the graft should be used for injection if possible.

Blood Access

Best results are achieved by:

1. Leaving the graft in place for approximately two weeks prior to use. Immediate use may increase the risk of hematoma formation.
2. Inserting the blood access needle at a 20° to 45° angle, with the bevel up until the graft is penetrated, and then advancing the needle parallel to the graft. **Routine flipping of needles is not recommended.**
3. **ROTATING THE CANNULATION SITES.** Repeated cannulation in the same area leads to formation of hematoma or pseudoaneurysm.
4. **NOT** cannulating within a dialysis needle's length of the proximal and distal anastomosis. See **PRECAUTION #6**.
5. Strict adherence to aseptic technique to minimize infection.
6. Applying moderate digital pressure to the cannulation site after needle withdrawal. This compression assists in hemostasis. (Note: There should always be a pulse or thrill of near equal intensity proximal and distal to the area of compression. Ensure that hemostasis is maintained with the least amount of pressure required.) See **WARNING #1**.

References

1. Health and safety aspects of 'Fluon' polytetrafluoroethylene, Technical Service Note F10, Sixth Edition, Imperial Chemical Industries PLC, 1989.
2. "Perigraft Seromas, Complicating Arterial Grafts," Robert M. Blumenberg, M.D., et al., Surgery, Vol. 97, No. 2, February 1985.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Referanslar

1. Health and safety aspects of 'Fluon' polytetrafluoroethylene, Technical Service Note F10, Sixth Edition, Imperial Chemical Industries PLC, 1989.
2. "Perigraft Seromas, Complicating Arterial," Robert M. Blumenberg, M.D., et al., Surgery, Vol. 97, No. 2, February 1985.

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STERILE | **EO**

Sterilized By Using Ethylene Oxide

Stérilisé à l'oxyde d'éthylène
Sterilisiert mit Ethylenoxid
Sterilizzato mediante ossido di etilene
Esterilizado usando óxido de etileno
Gesteriliseerd met behulp van ethyleenoxide
Esterilizado por óxido de etileno
Αποστειρωμένο με αιθυλενοξείδιο
Steriliseret ved etylenoxid
Steriliserad med etylenoxid
Steriloitu etyleenioksidilla
Steriliseret med etylenoksid
Product sterylizowany tlenkiem etylenu
Etilénoxiddal sterilizálva
Sterilizováno etylenoxidem
Etilen oksit ile sterilize edilmiştir

Manufactured By

Fabriqué par
Hersteller
Prodotto da
Fabricado por
Vervaardigd door
Fabricado por
Κατασκευάζεται από την
Fremstillet af
Tillverkad av
Valmistaja
Tilvirket av
Producent
Gyártó
Vyrábí
Uretici

Attention, See Instructions For Use

Attention, voir le mode d'emploi
Achtung! Siehe Gebrauchsanweisung
Attenzione, leggere le istruzioni per l'uso
Atención, consulte la instrucciones de uso
Let op, zie Gebruiksaanwijzing
Atenção, Consulte as Instruções de Utilização
Προσοχή, βλ. Οδηγίες Χρήσης
Bemærk! Se Brugervejledning
Obs! Se bruksanvisningen
Huomio, Lue käyttöohjeet
OBS! Se bruksanvisningen
Uwaga: Należy zapoznać się z instrukcją użytkowania
Figyelem, lásd a használati utasítást
Pozor, viz pokyny k použití
Dikkat, Kullanma Talimatına Bakınız

EC | **REP**

EEA Authorized Representative

Représentant autorisé EEA
EEA Berechtigter Vertreter
Rappresentante autorizzato EEA
Representante autorizado en el EEE
Gemachtigde binnen de EER
Representante autorizado no EEE
Εξουσιοδοτημένος αντιπρόσωπος στον EOX
EØF-autoriseret repræsentant
Auktoriserad representant inom EES
ETA:n valtuutettu edustaja
Autorisert EØS-representant
Autoryzowany przedstawiciel dla krajów Europejskiego Obszaru Gospodarczego
EGT meghatalmazott képviselő
Autorizovaný zástupce pro EHP
Avrupa Ekonomik Alanı Yetkili Temsilcisi

Contents

Contenu
Inhalt
Contenuto
Contentido
Inhoud
Conteúdo
Περιεχόμενα
Innhöld
Innehåll
Sisältö
Innhöld
Zawartość
Tartalom
Obsah
İçerikler

NON | **PYROGENIC**

Non-Pyrogenic

Apyrogène
Pyrogenfrei
Apirogeno
Apirógeno
Niet-pyrogeen
Apirogénico
Μη πυρετογόνο
Apyrogen
Pyrogenfri
Pyrogeeniton
Ikke-pyrogen
Apirogenny
Pirogénmentes
Apyrogenní
Pirojenik değildir

Do Not Use If Package Is Damaged Or Opened

Ne pas utiliser si l'emballage est ouvert ou endommagé
Nicht verwenden, wenn die Verpackung geöffnet oder beschädigt ist
Non utilizzare se la confezione è aperta o danneggiata
No usar si el envase está abierto o dañado
Niet gebruiken wanneer de verpakking geopend of beschadigd is
Não utilizar se a embalagem estiver aberta ou danificada
Μη χρησιμοποιείτε το προϊόν εάν η συσκευασία του έχει ανοιχθεί ή έχει υποστεί ζημιά
Må ikke bruges, hvis emballagen er beskadiget eller åbnet
Använd inte produkten om förpackningen öppnats eller skadats
Ei saa käyttää, jos pakkaus on auki tai vaurioitunut
Skal ikke brukes dersom pakken er skadet eller åpen
Nie stosować, jeśli opakowanie jest otwarte lub uszkodzone
Ne használja fel, ha a csomagolás nyitott vagy sérült
Nepoužívejte, pokud je obal poškozen či otevřen
Paket açılmış ya da hasarlıysa kullanmayınız

Variously protected by one or more of the following U.S. patents: 5,641,443; 5,827,327; 6,190,590; 6,203,735; 6,273,912; 6,436,135; 6,746,480. Other U.S. and foreign patents pending. Copyright ©2007 C. R. Bard, Inc. All rights reserved. Printed in the USA.

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Protection variable par un ou plusieurs des brevets américains suivants : 5.641.443 ; 5.827.327 ; 6.190.590 ; 6.203.735 ; 6.273.912 ; 6.436.135 ; 6.746.480. Autres brevets en instance aux États-Unis et dans d'autres pays. Copyright ©2007 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Geschützt durch ein oder mehrere der folgenden US-Patente: 5,641,443; 5,827,327; 6,190,590; 6,203,735; 6,273,912; 6,436,135; 6,746,480. Weitere Patente in den USA und anderen Ländern angemeldet. Copyright ©2007 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Protetto in varia misura da uno o più dei seguenti brevetti statunitensi: 5.641.443; 5.827.327; 6.190.590; 6.203.735; 6.273.912; 6.436.135; 6.746.480. Altre domande di brevetti statunitensi e stranieri in corso. Copyright ©2007 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Protegido por una o más de las siguientes patentes estadounidenses: 5.641.443; 5.827.327; 6.190.590; 6.203.735; 6.273.912; 6.436.135; 6.746.480. Otras patentes estadounidenses y extranjeras pendientes. Copyright ©2007, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Kunnen zijn beschermd door een of meer van de volgende Amerikaanse octrooien: 5,641,443; 5,827,327; 6,190,590; 6,203,735; 6,273,912; 6,436,135; 6,746,480. Andere Amerikaanse en buitenlandse octrooien aangevraagd. Copyright ©2007 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.
Protegido de diversas maneiras por uma ou mais das seguintes patentes americanas: 5,641,443; 5,827,327; 6,190,590; 6,203,735; 6,273,912; 6,436,135; 6,746,480. Outras patentes americanas e estrangeiras pendentes. Copyright ©2007, C. R. Bard, Inc. Todos os direitos reservados. Impresso nos E.U.A.