# *Instructions for Use*

***Venaflor® II ePTFE Vascular Graft  
Instructions For Use***

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**PERIPHERAL  
VASCULAR**

# INSTRUCTIONS FOR USE

**Caution:** Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

## Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

### Device Description

**Venaflo® II** ePTFE Vascular Grafts are constructed of expanded polytetrafluoroethylene (ePTFE) with carbon impregnated into the inner portions of the graft wall. This graft is also modified with a proprietary cuff with trim lines at the venous end.

Stepped configurations may help minimize the risk of arterial steal syndrome and high cardiac output. **CenterFlex™** graft configurations have an external spiral bead support on the center of the graft and can be used where resistance to compression or kinking is desired.

### Indications for Use

**Venaflo II** Vascular Grafts in **Straight, Stepped, CenterFlex,** and **Stepped CenterFlex** configurations are indicated for use as subcutaneous arteriovenous conduits for blood access only.

### Contraindications

None known.

### Warnings

1. All **Venaflo II** ePTFE vascular grafts are supplied sterile and non-pyrogenic unless the package is opened or damaged. **Venaflo II** grafts are sterilized by ethylene oxide. Each graft is intended for single patient use only. **DO NOT RESTERILIZE.**
2. Do not use after expiration date printed on the label.
3. Do NOT cannulate the externally supported portion of any **Venaflo II** ePTFE vascular graft due to the risk of beading embolization and/or pseudoaneurysm.
4. **DO NOT REMOVE THE EXTERNAL SPIRAL BEAD SUPPORT FROM ANY CenterFlex CONFIGURED GRAFT.** Attempts to remove the beading may damage the graft. If damage occurs, discard the graft.
5. **Venaflo II** grafts do not stretch (are non-elastic) in the longitudinal direction. Failure to cut the graft to an appropriate length may result in anastomotic or graft disruption, leading to excessive bleeding, loss of limb or limb function, and/or death.
6. Aggressive and/or excessive graft manipulation when tunneling, or placement within a too tight or too small tunnel, may lead to separation of the spiral beading and/or graft breakage. **DO NOT** pass the cuff portion (venous end) of the **Venaflo II** graft through a tunneler sheath or the tissue tunnel, as this could lead to separation of the spiral beading and/or graft breakage.
7. When embolectomy or balloon angioplasty catheters are used within the lumen of the graft, the inflated balloon size must match the inner diameter of the graft. Over-inflation of the balloon or use of an inappropriately sized balloon may dilate or damage the graft.
8. Avoid repeated or excessive clamping at the same location on the graft. If clamping is necessary, use only atraumatic or appropriate vascular smooth jawed clamps to avoid damage to the graft wall. Do not clamp the cuffed portion of the graft.
9. Exposure to solutions (e.g., alcohol, oil, aqueous solutions, etc.) may result in loss of the graft's hydrophobic properties. Loss of the hydrophobic barrier may result in graft wall leakage. Preclotting of this graft is unnecessary.<sup>1</sup>
10. Avoid excessive graft manipulation after exposure to blood or body fluids. Do not forcibly inject any solution through the lumen of the graft, or fill the graft with fluid prior to pulling it through the tunnel as loss of the graft's hydrophobic properties may occur. Loss of the hydrophobic barrier may result in graft wall leakage.
11. Do NOT expose **Venaflo II** grafts to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.<sup>2</sup>
12. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.
13. During tunneling, be sure to create a tunnel that closely approximates the outer diameter of the graft. A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation.
14. Care should be taken when cannulating the graft for dialysis access; avoid the cuff area and the **CenterFlex** external support, properly rotate cannulation sites, and use proper compression to achieve hemostasis during post-cannulation care.

### Precautions

1. Only physicians qualified in vascular surgery techniques should implant this prosthesis. The healthcare provider is responsible for all appropriate postoperative care instructions to the patient.

2. The healthcare provider must observe aseptic technique during implantation, postoperatively, and during cannulation.
3. When suturing, avoid excessive tension on the suture line, inappropriate suture spacing and bites, and gaps between the graft and host vessel. Failure to follow correct suturing techniques may result in suture-hole elongation, suture pull-out, anastomotic bleeding and/or disruption. Refer to "Suturing" for further instructions.
4. Do not cannulate the cuffed portion of the graft. For further instructions, refer to the "Blood Access" section.
5. Consider intraoperative and postoperative anticoagulation therapy for each patient as appropriate.

### Adverse Reactions

Potential complications that may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: Disruption or tearing of the suture line, graft, and/or host vessel; suture-hole bleeding; graft redundancy; thrombosis; embolic events; occlusion or stenosis; ultrafiltration; seroma formation; swelling of the implanted limb; formation of hematomas or pseudoaneurysm; infection; arterial steal syndrome; aneurysm/dilation; blood leakage; hemorrhage; and/or skin erosion.

# DIRECTIONS FOR USE

### Equipment Required

Tunneler, double armed suture, atraumatic clamp, scalpel and/or scissors.

### Opening the Package

Hold the outer tray in one hand. Peel back the lid. Remove the inner tray. Peel back the inner tray lid slowly and carefully remove the graft using sterile atraumatic instruments or gloves. Protect the graft against damage from sharp or heavy instruments.

# General Operative Techniques

### Sizing



Use the following guidelines when sizing the cuffed end of all **Venaflo II** ePTFE vascular grafts for the venous anastomosis.

Size of Vein (External Diameter)	Cuff Sizing (See Figure 1)
>8mm	Do not trim
6-8mm	Trim on the line labeled 6-8
3-5mm	Trim on the line labeled 3-5

### Trimming Recommendations

Hold the graft with the cuff in an open configuration, as shown in **Figure 1**.

Using curved scissors, start trimming the cuff on the dashed line located nearest the heel. (See **Figure 2**.) Continue trimming on the dashed lines following the curvature suggested. **Note: Do not fold the graft cuff in half when trimming. Do not trim material from the heel.**

### Tunneling Recommendations

Create a tunnel that closely approximates the outer diameter of the graft. Always follow the Instructions for Use for the specific tunneler utilized to place the **Venaflo II** ePTFE vascular graft.

**WARNING:** A tunnel that is too loose may result in delayed healing and may also lead to perigraft seroma formation. See **Warning #13**.

**Note:** Prior to utilizing a sheath tunneler, verify that the **Venaflo II** graft O.D. dimensions fit the sheath I.D. dimensions.

### Venaflo II ePTFE Vascular Grafts in Straight Anatomical Configurations:

After selecting the appropriate tunneler, introduce the tunneler through the arterial site. After creating the tunnel, feed the graft from the venous end to the arterial end. (See **Figure 3**.) Ensure the cuffed end does not pass through the tunnel. Only the arterial end of the **Venaflo II** ePTFE vascular graft should be sized for length.

Figure 3

**Venaflo® II ePTFE Vascular Grafts in Looped Anatomical Configurations:**

Perform an incision to expose the artery and vein to be anastomosed and a counter-incision at the site where the graft will curve. For **CenterFlex™** grafts, this is where the flexed portion of the graft should be placed. Position the graft overlying the forearm to optimize graft placement prior to tunneling. Introduce the tunneler at the counter-incision site.

After creating the first tunnel, replace the bullet tip (if necessary), attach the arterial end of the graft to the tunneler, and pull the graft through the tunnel leaving the cuffed end exposed near the vein to be anastomosed.

For the second tunnel, repeat the tunneling procedure, initiating the tunnel at the anastomotic site and advancing it down to the counter incision. (See Figure 4.) Once again, ensure that the venous (cuffed) graft end does not pass into the tunnel.

Figure 4

Only the arterial end of the **Venaflo II** ePTFE vascular graft should be sized for length.

**Suturing**

Size the graft appropriately to minimize excessive tension at the suture line. Use a tapered, non-cutting needle with a non-absorbable monofilament suture approximately the same size as the needle. Take 2 mm suture bites in the graft following the curve of the needle and gently pull the suture at a 90° angle. Proper sizing of the graft length prior to implant will minimize suture-hole elongation caused by excessive tension. See **Precaution #3**.

**Suggested Suturing Technique**

After sizing the cuff to the vein diameter, perform a venotomy five to eight millimeters shorter than the length of the cuff.

Using double-armed suture, tack the heel of the cuff, aligning the toe in the direction of venous blood flow.

Suture continuously from the heel to just beyond the midpoint of the cuff on one side of the venotomy. Using the other arm of the heel suture, continue suturing to just beyond the midpoint of the graft on the other side of the venotomy.

Extend the venotomy to the exact cuff length. (See Figure 5.)

Figure 5

Complete the "cuffed" anastomosis by suturing with one of the suture arms, around the toe, to the other side of the cuff. Always tie off on the side of the cuff (See Figure 6).

Figure 6

Perform the arterial anastomosis using standard surgical techniques. In order to retain a cuffed configuration, only the arterial end of the **Venaflo II** ePTFE vascular graft should be sized for length.

**Thrombectomy**

Techniques for declotting **Venaflo II** ePTFE vascular grafts include but are not limited to the use of balloon catheters.

In the event of a postoperative occlusion, **Venaflo II** grafts can be declotted (see Warning #7) using either a longitudinal or transverse incision. During the early postoperative period, the natural progression of wound healing renders the graft translucent in appearance. When the material is in this state, it is recommended that a longitudinal incision with stay sutures be used for thrombectomy procedures. If a transverse incision is performed, the use of a horizontal mattress suture technique may be required to aid in closure.

**Longitudinal Incision:**

Place stay sutures before introducing the embolectomy catheter. (If the graft is a **CenterFlex** product, cut through the spiral beading and base tube. After closure, the spiral beading will realign itself.) Place a longitudinal incision in the graft that is long enough to accommodate the extraction of a fully dilated thrombectomy catheter balloon. A patch may be considered as an aid to graft closure.

**Transverse Incision:**

No stay sutures are necessary. A horizontal mattress suture is recommended for graft closure.

**Angiography**

Should angiography be performed at the time of the procedure, the artery proximal to the graft should be used for injection, if possible.

**Blood Access**

Best results are achieved by:

1. Leaving the graft in place for approximately two weeks prior to use. Immediate use may increase the risk of hematoma formation.
2. Inserting blood access needle at a 20° to 45° angle, with the bevel up until the graft is penetrated, and then advancing the needle parallel to the graft. Routine flipping of needles is not recommended.
3. ROTATING THE CANNULATION SITES. Repeated cannulation in the same area may lead to formation of hematoma or pseudoaneurysm.
4. Cannulating greater than a dialysis needle's length of the proximal anastomosis. See **Precaution #4**.
5. Strictly adhering to aseptic technique to minimize infection.
6. Applying moderate digital pressure to the cannulation site after needle withdrawal. This compression aids in hemostasis.

**Note: There should always be a pulse or thrill near equal intensity proximal and distal to the area of compression. Ensure that hemostasis is maintained with the least amount of pressure required.**

**References**

1. "Perigraft Seromas, Complicating Arterial Grafts," Robert M. Blumenberg, M.D., et al., Surgery, Vol. 97, No. 2, February 1985.
2. *Guide to the Safe Handling of Fluoropolymer Resins*, 3rd Edition, The Fluoropolymers Division of the Society of the Plastics Industry, Inc.

**Warranty**

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Variouly protected by one or more of the following U.S. Patents:****5,827,327, 6,190,590, 6,203,735, 6,273,912, 6,436,135 and 6,746,480.****Other U.S. and Foreign Patents Pending.****Copyright © 2008 C. R. Bard, Inc. All Rights Reserved. Printed in the U.S.A.**
Copyright © 2008 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2008 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2008 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2008 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright © 2008 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.