Instructions For Use

Bard® DeBakey® Soft Woven Polyester Vascular Prostheses. 1

Mode d'emploi

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Bard® DeBakey® Soft Woven Polyester Vascular Prostheses

Caution: Federal (U.S.A.) law restricts these devices to sale, distribution and use by or on the order of a physician. Only physicians qualified in vascular surgery techniques should use these prostheses.

Instructions For Use

Product Description:

These vascular prostheses are manufactured with PET (Polyethylene terephthalate) polyester fibers and feature crimped construction to facilitate graft placement.

Indications for Use:

These woven vascular prostheses are indicated for replacement or bypass procedures in aneurysmal and occlusive disease of the thoracic and abdominal aorta, visceral arteries and proximal peripheral arteries, exclusive of the coronary arteries. Straight woven prostheses are also indicated for use in intra-aortic balloon assist procedures.

Contraindications:

These prostheses are contraindicated for use in coronary vascular procedures and for use as an arteriovenous conduit for blood access.

Warnings:

1. **These prostheses must be properly preclotted with non-heparinized blood prior to exposure to full arterial pressure to avoid unnecessary bleeding or blood loss.**
2. **While it is common practice to employ moderate doses of intraoperative heparin, care should be taken not to exceed the manufacturer's recommended dose for such procedures. Excessive amounts can result in bleeding.**
3. **Due to variability in patient response to heparinization, it is essential that the adequacy of anticoagulation during surgery and the precision of neutralization at the conclusion of the procedure be closely monitored. Adherence to a strict protocol, determined by each hospital, can prevent excessive bleeding.**
4. **As with any vascular prosthesis, occasional difficulties with hemostasis may occur. In the event that hemostasis is not easily obtained, the surgeon may wish to consider the following:**
 - Investigate for systemic coagulopathy and treat appropriately.
 - Compression, as necessary.
 - Additional sutures and/or pledgets at the anastomosis, as necessary.
 - Reversal of heparin with protamine sulfate, as necessary.
 - Utilization of topical coagulation therapy, e.g., thrombin, as necessary.

Precautions:

1. STERILE, unless the package is opened or damaged. SINGLE USE ONLY.
2. Bard DeBakey woven vascular prostheses are sterilized by ethylene oxide. Resterilization is not recommended.
3. Care should be taken when clamping the prosthesis to avoid damage to the fibers and preclot. This risk may be reduced with the use of atraumatic, shod-type vascular clamps.
4. Taper point suture needles are recommended so that the prosthesis fibers will not be cut.
5. To reduce the risk of fraying, care should be taken when cutting woven prostheses. It is recommended that a hot cautery be used.
6. To minimize fluid collection around the graft in aorto-femoral and femoral popliteal reconstructive procedures, the lymphatics should be carefully ligated and sealed, especially in the groin area.
7. In the abdomen, care should be taken to separate the prosthesis and anastomoses from adjacent bowel structures by protecting with an adequate layer of retroperitoneal tissue, or omentum / peritoneum to reduce the risk of secondary fistula formation or bowel erosion.
8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Adverse Reactions:

Adverse reactions that may occur with the use of these products or with any vascular graft include perioperative hemorrhage, thrombosis, thromboembolic complications, infection, anastomotic aneurysms, and aortoenteric fistula formation.

Equipment Required:

Cautery
Non Resorbable suture with taper point needles
Atraumatic shod-type clamps

Directions for Use

Graft Preparation and Preclotting:

Even though these prostheses are manufactured to have low permeability, preclotting is recommended. Precise preclotting is necessary if there is evidence of coagulopathy, or in the presence of high doses of heparin.

A suggested method of preclotting is described below. However, it is appropriate to employ whatever method your institution has accepted and validated as standard practice for woven grafts.

1. Withdraw non-heparinized blood (25-35 ml) from the patient. Thoroughly wet the graft with blood. Clamp the distal openings of the graft and force blood through the interstices of the graft until the interstices are sealed and the graft is adequately preclotted. A second syringe of blood may be necessary. Allow the prosthesis to sit in the blood for several minutes.
2. With fingers, strip the excess blood and clot from the prosthesis and gently suction the lumen. The prosthesis is now ready for use.
3. After the proximal anastomosis of the prosthesis has been completed, it is suggested that the graft preclotting be tested while checking for anastomotic leaks. Occlude the distal end of the graft while releasing the occluding clamp on the artery. Observe the graft and anastomosis. If bleeding is noted, see Warnings for additional hemostasis information.

WARRANTY

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issued or revision date and a revision number for these instructions are included for the user's information. In the event that 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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**Bard and DeBakey are registered trademarks of C. R. Bard, Inc. or an affiliate.**
**Copyright © 2008 C. R. Bard, Inc. All rights reserved. Printed in the USA.**
Copyright © 2008 C. R. Bard, Inc. Tous droits r serv s. Imprim  aux  tats-Unis.
Copyright © 2008 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2008 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2008, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.