
# Instructions for Use

## **Bard® Parsonnet™ Vascular Probes**

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# Bard® Parsonnet™ Vascular Probes

## Instructions For Use

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

### Product Description:

The **Parsonnet** Vascular Probe is a sterile, single-use device. The shaft of the probe is constructed with braided textiles, covered with a smooth, radiopaque polyurethane coating. The shafts have printed graduations at 1cm intervals to facilitate localization of obstructions. The probe tip diameters are available in .5mm increments, from 1.0mm to 3.5mm. Each probe is equipped with a removable, flexible, stainless steel stylet, which provides extra stiffness to the shaft of the probe, if necessary when intubating an artery.

### Indications for Use:

**Parsonnet** Vascular Probes are designed for use during coronary artery and peripheral vascular surgery. They are used to assist in suture placement and maintenance of a dry operative field by preventing antegrade and retrograde bleeding during anastomosis. They are used for calibration of the internal diameter of the artery, and for location and calibration of stenotic lesions remote from the anastomotic site.

### Contraindications:

None known.

### Warnings:

1. **To reduce the risk of vessel wall damage, care should be taken when inserting the probe into a vessel. Choose the correct size probe for the vessel. Do NOT force the tip into the lumen of the vessel. If resistance is encountered on intubation, remove the probe. Determine the cause of the resistance and take appropriate action, such as using a smaller size probe.**
2. **When an endarterectomy has been performed, the probes must be inserted with care, to reduce the risk of inadvertent intimal dissection.**

### Precautions:

1. This device should be used only by physicians qualified in cardiovascular and peripheral vascular surgery techniques.
2. STERILE, unless the package is opened or damaged. Single use only.
3. These products are sterilized by ethylene oxide. DO NOT RESTERILIZE.
4. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

### Adverse Reactions:

Adverse reactions that may occur with the use of this product or with any vascular probe include vessel wall damage, intimal dissection, or disruption of plaque, with resultant embolization.

## Directions For Use:

### Coronary Artery Surgery

1. The diseased coronary artery to be bypassed is opened longitudinally in the area chosen to receive the bypass graft. Once the arteriotomy has been made, probes are inserted proximally and distally into the recipient artery, either one at a time or simultaneously. The probes are removed after both corners and one side of the anastomosis have been completed **and after** reestablishment of coronary flow (if ischemic arrest has been used). After the removal of the probe, the lumen is kept open with microforceps to allow free escape of any entrapped air. Once air is removed and flow is established, final sutures are placed across the bleeding edges to complete the anastomosis.
2. If ischemic arrest is not used, the probes may be used to reduce bleeding into the anastomotic site by being placed within the vessel in an antegrade and/or retrograde fashion.

### Peripheral Vascular Surgery

1. After an arteriotomy has been made, the appropriately sized probe may be inserted proximally and/or distally to determine vessel size, patency and location of stenoses.
2. The probe may be used to determine the course of the vessel and for location of major tributaries.

### Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

### Reference:

Parsonnet, V., Gilbert, L., Gielchinsky, I.: Graduated probes for coronary bypass surgery. *The Journal of Thoracic and Cardiovascular Surgery*, 68(3):424-427, 1974.

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#### Bard and Parsonnet are Trademarks and/or Registered Trademarks of C. R. Bard, Inc. or an affiliate
**Copyright ©2008 C. R. Bard, Inc. All rights reserved. Printed in the USA.**
Copyright ©2008 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright ©2008 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright ©2008 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright ©2008 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.