RIVAL™

PTA Balloon Dilatation Catheter

ENGLISH

INSTRUCTIONS FOR USE

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.

Device Description

The RIVAL™ PTA Balloon Dilatation Catheter is a balloon catheter consisting of an over the wire catheter with a balloon fixed at the distal tip. The non-compliant nylon, low profile balloon is designed to provide consistent balloon diameters and lengths. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. The proprietary balloon technology and .035" guidewire compatible dual lumen catheter are designed to optimize the balance between pushability and trackability. The proximal portion of the catheter includes an inflation female luer lock hub and a guidewire female luer-lock hub. Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use. A re-wrapping tool is also provided on the catheter shaft. A stylet is placed into the tip of the catheter to aid in rewrap/refolding of the balloon. This product is not manufactured with any latex.

Indications for Use

The RIVAL™ PTA Balloon Dilatation Catheter is intended to dilate stenosis in the peripheral arteries, treat obstructive lesion of native or synthetic A-V fistulae, and/or re-expand endoluminal stent graft elements in the iliac arteries.

Contraindications

None known.

Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Do not reuse, reprocess or re-sterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not re-sterilize. After re-sterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing and/or re-sterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
6. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent overpressurization, use of a pressure monitoring device is recommended.
7. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The RIVAL™ PTA Balloon Dilatation Catheter shall only be used by physicians trained in the performance of percutaneous transluminal angioplasty.
3. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size sheath introducer than indicated on the label.
4. Use the recommended balloon inflation medium (30-50% contrast medium/50-70% sterile saline solution). It has been shown that a 30/70% contrast / saline ratio has yielded faster balloon inflation / deflation times. Never use air or other gaseous medium to inflate the balloon.
5. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
6. If resistance is felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter, introducer sheath and guidewire (as necessary) as a single unit.
7. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
8. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.

Potential Adverse Reactions

The complications which may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

Directions for Use

Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources. Rotate inventory so that the catheters and other dated products are used prior to the 'Use By' date. Do not use if packaging is damaged or opened.

Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- 0.035" guidewire

Dilatation Catheter Preparation

1. Remove Catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the balloon guard and stylet by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon protector and stylet and slide distally off of the balloon catheter.
3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (30-50% contrast medium/50-70% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 as needed until bubbles no longer

appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.

8. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

Use of the RIVAL™ PTA Balloon Dilatation Catheter

1. Backload the distal tip of the RIVAL™ PTA balloon catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and inflate the balloon to the appropriate pressure.
4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy and that no contrast is left in the balloon.
5. While maintaining negative pressure and the position of the guidewire, grasp the balloon catheter just outside the sheath and withdraw the deflated dilatation catheter over the wire through the introducer sheath. Use of a gentle counterclockwise motion may be used to help facilitate catheter removal through the introducer sheath.

Balloon Reinsertion

Precaution: Do not continue to use the balloon catheter if the shaft has been bent or kinked.

Precaution: Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire or stylet.

1. Load the balloon catheter onto a guidewire or stylet.
2. Advance the balloon re-wrap tool over the catheter to the proximal end of the balloon.
3. Grasp the catheter shaft just proximal to the balloon with one hand, and with the other hand gently slide the rewrap tool over the balloon to the catheter tip and then back over the balloon to the catheter.
4. Slide the re-wrap tool to the proximal end of the catheter shaft.
5. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
6. Continue the procedure according to the "Use of the RIVAL™ PTA Balloon Dilatation Catheter" section herein.

Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Bard and Rival are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.
Copyright © 2009 C. R. Bard, Inc. All rights reserved.
Copyright © 2009 C. R. Bard, Inc. Tous droits réservés.
Copyright © 2009 C. R. Bard, Inc. Alle Rechte vorbehalten.
Copyright © 2009 C. R. Bard, Inc. Tutti i diritti riservati.
Copyright © 2009 C. R. Bard, Inc. Reservados todos los derechos.