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IMPRA[®]
Innovation for Optimal Outcomes[™]

Peripheral Vascular Sheath Tuner and Sterilization Cassette

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

INFORMATION FOR USE


Indication

The Peripheral Vascular Sheath Tuner (PV Sheath Tuner) is intended to aid in creating subcutaneous tunnels for the placement of vascular prostheses (grafts) for peripheral vascular, femoral arteriovenous access, and extra-anatomic bypass procedures. The tuner will facilitate the placement of externally supported grafts up to 8mm I.D. A Tuner Sterilization Cassette is provided to hold and protect the Tuner components when not in use and to allow penetration of steam during sterilization. The Cassette is to be used with an appropriate sterilization wrap.

Description

The PV Sheath Tuner consists of four components, that when assembled create a single unit device. The Tuner components include the Handle, Inner Rod, Outer Sheath, and Bullet Tip (reference Figure 1).

Figure 1

Both the Tuner and the Tuner Sterilization Cassette are reusable medical devices. The Tuner components are made of medical grade stainless steel. The Tuner Sterilization Cassette is composed of Radel® R (polyphenylsulfone) and has four medical grade stainless steel latches. These materials are designed to withstand repeated use, cleaning, and resterilization.

Precautions

1. Prior to utilizing the PV Sheath Tuner, verify that the vascular graft outside diameter fits within the Sheath and the graft length is longer than the tissue tunnel to be created. Failure to do so may make it difficult to achieve optimal graft placement.
2. Care should be exercised to protect the Tuner from mechanical damage. Rough handling may mar the surface of the Tuner.
3. Do not soak Tuner components in saline or bleach solutions. Doing so could cause damage to component surfaces.
4. The Tuner is supplied non-sterile and must be thoroughly cleaned and sterilized prior to each use. Cleaning/sterilization of the Tuner prior to use is the responsibility of the healthcare facility.

Cleaning Prior To Sterilization

Personnel should wear appropriate protective attire when handling the PV Sheath Tuner components and the Tuner Sterilization Cassette. Immediately after use, completely disassemble the Tuner and rinse the components and the Tuner Sterilization Cassette in cool running tap water, utilizing brushing if needed, to aid in the removal of visible debris. Use stainless steel wire or stiff plastic brushes that are compatible with stainless steel surgical instruments. A mechanical cleaning system or manual cleaning method may be used to clean the Tuner components and the Tuner Sterilization Cassette. Consult the mechanical cleaning system and cleaning agent manufacturer's instructions for cleaning surgical instruments.

The following manual cleaning method is suggested:

1. Remove all particles of adhering tissue and dried blood from the Tuner, Sterilization Cassette and the Tuner components using cool running tap water and mild, enzymatic detergent. Prepare the enzymatic detergent and soak the Tuner components and Cassette per the detergent manufacturer's instructions. Use stainless steel wire or stiff plastic brushes that are compatible with stainless steel surgical instruments to aid in cleaning the inner lumen of the Sheath and crevices on the threaded portions of the Rod, Bullet Tip, and Handle.

2. Immediately after washing the Tunneler components and the Tunneler Sterilization Cassette, rinse in distilled water to remove all traces of detergents and debris. Allow all Tunneler components and the Tunneler Sterilization Cassette to air dry.
3. It is important to inspect all Tunneler components and the Tunneler Sterilization Cassette to insure that they are clean and free of damage prior to sterilization. Damage that could prevent proper function of the Tunneler includes scratches, cross-threading, etching, pitting, burrs, nicks, dents or cracks. For proper use, the Tunneler Sterilization Cassette should latch when closed and should be free of any structural damage. **Reference PRECAUTIONS #2, #3, and #4.**

Sterilization

After cleaning and inspection, all disassembled PV Sheath Tunneler components should be placed in the Tunneler Sterilization Cassette. Prior to use, the Tunneler Sterilization Cassette must be enclosed in an appropriate commercially available sterilization wrap, used in accordance with the wrap manufacturer's instructions.

Warnings:

1. **Do not stack Tunneler Sterilization Cassettes during sterilization period.**
2. **Do not sterilize the Tunneler in the Tunneler Sterilization Cassette using a 250°F (121°C) gravity displacement cycle.**

Sterilize the Tunneler using the following cycle parameters:

Using a gravity displacement steam sterilizer:

Temperature: at or above 270°F (132°C)

Exposure Time: 15 minutes

Drying Time: 45 minutes

Using a pre-vacuum (high vacuum) steam sterilizer:

Temperature: at or above 270°F (132°C)

Exposure Time: 3 minutes

Drying Time: 45 minutes

Prior to opening the sterilization wrap, inspect the wrap for tears, pinholes, or other damage that might compromise sterility. The Cassette must be enclosed in an appropriate commercially available sterilization wrap in order to preserve sterility after processing. The following standards/guideline may be referenced for further information on the safe handling and biological cleaning and sterilization of medical devices:

"Recommended Practices for Cleaning and Caring for Surgical Instruments and Powered Equipment," AORN Journal, March 2002, Vol. 75, No. 3, pp. 627-641.

AAMI/ANSI/ISO 11134 "Sterilization of health care products - Requirements for validation and routine control - Industrial moist heat sterilization."

ANSI/AAMI ST35-1996 "Good Hospital Practice: handling and biological decontamination of safe medical devices in health care facilities and in non-clinical settings."

Using the Tunneler

1. Remove the sterile PV Sheath Tunneler components from the Tunneler Sterilization Cassette using appropriate aseptic technique.
2. To assemble the sterile Tunneler components: Carefully secure (thread) the Handle onto either end of the Rod. Slide the Rod through the Sheath, so that the threaded end of the Rod extends beyond the Sheath end. Attach (thread) the Bullet Tip to the exposed end of the Rod. Hand-tighten the assembly, making sure that the shoulders of the Bullet Tip and Handle are snug against the ends of the Sheath.
3. Follow standard surgical practice to create entrance and exit incisions for the Tunneler.
4. Use the sterile assembled Tunneler to create a tissue tunnel connecting the entrance and exit incisions.
5. When the Tunneler passes through the exit incision site and a satisfactory tunnel has been created, unscrew the Handle 1/2 turn, and then completely remove the Bullet Tip from the Rod.

- When tunneling is completed and the Bullet Tip will not be needed for the remainder of the procedure, remove the Bullet Tip from the surgical field and return it to the Tunneler Sterilization Cassette.
- Secure (tie) the graft onto the Rod by placing a suture through the tapered hole at the exposed end of the Rod. Verify that the vascular graft outside diameter fits within the Sheath and the graft length is longer than the tissue tunnel.
Reference PRECAUTION #1.
- Carefully pull the graft through the Sheath, using the Handle and Rod assembly.
 - When the graft has been pulled into position within the Sheath, release it from the Rod by cutting the suture.
 - While holding the graft at the exit incision site, grasp the protruding Sheath end and carefully pull (remove) the Sheath from the tissue tunnel, leaving the graft properly placed inside the tissue tunnel.
- Proceed with the vascular anastomoses.

Maintenance Procedures

Proper maintenance of the PV Sheath Tunneler requires careful handling and inspection of each Tunneler component and the Tunneler Sterilization Cassette. Particular care should be taken to avoid damaging the threaded portions of the Rod, Handle, and Bullet Tip. If these areas are damaged through improper handling (i.e., dropping, cross-threading, brushing and/or cleaning with abrasives) the parts may not fit together properly or be inter-changeable on the Rod ends. Do not use the Tunneler if such damage occurs. Contact Bard Peripheral Vascular for appropriate replacement parts and ordering information.

Warranty

Bard Peripheral Vascular, Inc. warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states or countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state or country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Bard and **IMPRA** are registered trademarks of C. R. Bard, Inc., or an affiliate.

Innovation for Optimal Outcomes is a trademark of C. R. Bard, Inc., or an affiliate.

Radel is a registered trademark of Solvay Advanced Polymers, LLC.

Copyright ©2004, C. R. Bard, Inc. All rights reserved. Printed in U.S.A.

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"Recommended Practices for Cleaning and Caring for Surgical Instruments and Powered Equipment," AORN Journal, maart 2002, vol. 75, nr. 3, p. 627 - 641.

AAAMI/ANSI/ISO 11134 "Sterilization of health care products - Requirements for validation and routine control – Industrial moist heat sterilization."

ANSI/AAAMI ST35-1996 "Good Hospital Practice: handling and biological decontamination of safe medical devices in health care facilities and in non-clinical settings."

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Copyright ©2004, C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright ©2004, C. R. Bard, Inc.. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright ©2004, C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti
Copyright ©2004, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.
Copyright ©2004, C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.
Copyright ©2004, C. R. Bard, Inc. Todos os direitos reservados. Impresso nos E.U.A.