ePTFE Cardiovascular Patch

Instructions for Use

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
PERIPHERAL
VASCULAR

INSTRUCTIONS FOR USE

Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

Device Description, Indications, Contraindications, Warnings, Precautions and Adverse Reactions

Device Description

The Expanded PTFE (ePTFE) Cardiovascular Patch is constructed from expanded polytetrafluoroethylene. All Cardiovascular Patches are supplied sterile and non-pyrogenic unless the package is opened or damaged. This is a single use device.

Indications for Use

The Cardiovascular Patch is indicated for use in repair and closure of the cardiovascular system.

Contraindications

The Cardiovascular Patch is NOT indicated for:

- Reconstruction of hernias and tissue deficiencies.
- Reconstruction and repair of passive biological membranes such as dura mater, pericardium, or peritoneum.

Use of this prosthesis in noncardiovascular applications can cause potentially serious complications, such as suture pull out, failure of closure or repair, or undesired healing to surrounding tissues.

Warnings

1. **Cardiovascular Patches are supplied sterile and non-pyrogenic unless the package is opened or damaged. Cardiovascular Patches are sterilized by ethylene oxide. Each patch is intended for single patient use only. DO NOT RESTERILIZE.**
2. **Do not use after expiration date printed on the label.**
3. **Cardiovascular Patches are not elastic. Care must be taken to cut the patches sufficiently large enough to eliminate anastomotic stresses.**
4. **Exposure to solutions (e.g., alcohol, oil, aqueous solutions, etc.) may result in loss of the patch's hydrophobic properties. Loss of the hydrophobic barrier may result in patch wall leakage. Preclotting of this patch is unnecessary.**
5. **Avoid excessive patch manipulation or pressure after exposure to blood or body fluids, as loss of the patch's hydrophobic properties may occur. Loss of the hydrophobic barrier may result in patch leakage.**
6. **After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.**

Precautions

1. Only physicians qualified in cardiovascular surgery techniques should use this prosthesis. The healthcare provider is responsible for all appropriate postoperative care instructions to the patient.
2. When suturing, avoid excessive tension on the suture line, inappropriate suture spacing and bites, and gaps between the patch and host vessel. Failure to follow correct suturing techniques may result in suture hole elongation, suture pull-out, anastomotic bleeding and/or disruption. Refer to "Suturing" for further instructions.
3. DO NOT use absorbable sutures.
4. DO NOT use less than a 2 mm suture bite in cardiac and great vessel applications.
5. DO NOT use less than a 1 mm suture bite in peripheral vascular applications.
6. DO NOT suture two patches together to make a larger patch.
7. DO NOT contaminate or damage the patch.
8. DO NOT patch in areas that may be subjected to repeated punctures.

Adverse Reactions

Possible complications with the use of any cardiovascular patch include, but are not limited to: infection, thrombus formation, anastomotic blood leakage, peripatch seroma formation, pseudoaneurysm formation, and hematomas.

DIRECTIONS FOR USE

Equipment Required

Vascular forceps, needle holder, taper point, non-resorbable sutures, and topical hemostatic agent.

Opening the Package

Hold the outer tray in one hand. Peel back the lid. Remove the inner tray. Peel back the inner tray lid slowly and carefully remove the patch using sterile atraumatic instruments or gloves. Protect the patch against damage from sharp or heavy instruments. The patch is ready to be sized for repair or closure.

Sizing

It is critical that the patch be cut to the proper size for closure or repair. Cutting the patch too small increases the tension on the suture line which may result in elongated suture holes, thus causing increased suture line bleeding. Never suture two patches together to form a larger patch. This is not an accepted design criteria of the patch and may result in an inadequate repair or closure.

Suturing

Best results are achieved using a tapered, non-cutting needle with a non-absorbable monofilament suture approximately the same size as the needle. Follow the curve of the needle to minimize suture hole elongation and bleeding from suture holes. Application of a topical hemostatic agent may aid in controlling any bleeding that may occur. The manufacturer's instructions for these products should be followed.

When placing a suture, avoid excessive tension on the suture line and incorporate sufficient material in the stitch. Use a minimum 2 mm suture bite for cardiac and great vessel repairs; a minimum 1 mm suture bite for peripheral vascular applications.

Warranty

Bard Peripheral Vascular, Inc. warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular, Inc.'s sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and a revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular, Inc. to see if additional product information is available.

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**Bard and IMPRA are registered trademarks of C. R. Bard, Inc. or an affiliate.**
Copyright © 2008 C. R. Bard, Inc. All rights reserved. Printed in the U.S.A.
Copyright © 2008 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2008 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2008 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2008 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.