Instructions For Use	
Bard® Carotid Bypass Shunts	1
Mode d'emploi	
Shunts carotidiens Bard®	3
Gebrauchsanweisung	
Bard® Karotis-Bypass-Shunts	5
Istruzioni per l'Uso	
Shunt Bard® per Bypass Carotideo	7
Instrucciones de uso	
Shunts de bypass carotídeo Bard®	9
Gebruiksaanwijzing	
Bard® Carotis-bypass-shunts	11
Instruções de utilização	
Shunts de bypass carotídeo Bard®	13
Οδηγίες χρήσης	
Shunts Παράκαμψης Καρωτίδας Bard®	15
Brugervejledning	
Bard® Carotis Bypass Shunts	17
Bruksanvisning	
Bard® Carotis bypassshuntar	19
Käyttöohjeet	
Bard® -kaulavaltimon ohitussuntit	21
Bruksanvisning	
Bard® carotis bypass-shunter	23
Instrukcja użycia	
Połączenie pomostowa tętnicy szyjnej Bard®	25
Használati tájékoztató	
Bard® carotis bypass sönt	27
Pokyny k použití	
Shunty Bard® pro přemostění karotidy	29
Kullanma Talimatı	
Bard® Karotis Bypass Şantı	31

Bard® Carotid Bypass Shunts

Instructions For Use

Caution: Federal (U.S.A.) law restricts these devices to sale by or on the order of a physician. Only physicians qualified in vascular surgery techniques should use these products.

Product Description:

Bard Carotid Bypass shunts are cardiovascular cannula constructed of polyvinyl chloride. These shunts are sterilized by ethylene oxide and are sterile unless the package is opened or damaged.

- **Javid™** Carotid Shunts are tapered and long enough to allow for an external loop.
- Straight Carotid Shunts are available with or without beveled tips and side holes.
- **Brener™** Carotid Shunts are tapered with a side arm for aspiration.
- **Burbank™** Carotid Shunts are tapered, non-looping, and contain depth markings to facilitate placement.

Indications For Use:

Bard Carotid Shunts are designed to be a temporary blood conduit in the carotid arteries during a carotid endarterectomy procedure.

Contraindications:

Do not use a carotid bypass shunt if the carotid arteries demonstrate such diffuse atherosclerosis that it would prevent proper insertion and placement of the shunt.

Warnings:

1. **Single Use Only. DO NOT RESTERILIZE.**
2. **Assure that the shunts are properly stabilized in the artery or slippage may occur. Carefully read the Directions For Use.**
3. **These shunts are for temporary use and should not be permanently implanted.**
4. **Do not force a shunt that is too large into an artery. This may result in vessel disruption or damage.**

Precautions:

1. While blood is flowing through the shunt, clamp the side arms or T-extensions, unless flushing or aspirating the shunt.
2. Each shunt should be inspected prior to use to make certain that the lumen and eyes are clear and patent and the shunt has not been kinked or damaged during shipment or storage.
3. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state, and federal laws and regulations.

Adverse Reactions:

Adverse reactions that may occur with the use of these products or with any shunt used in a carotid endarterectomy procedure include: intimal or vessel wall disruption or damage, plaque or clot embolization, thrombosis, bleeding, or occlusion.

Equipment Required:

Fixation/Stabilization devices
Vessel loops
Umbilical-type tapes
Heavy ligatures
Tourniquets, or surgical clamps

Directions For Use

How Supplied/Stored:

Bard Carotid Bypass Shunts are supplied sterile, unless package is damaged or opened. Single use. Do not resterilize. **Bard** Carotid Bypass Shunts are single-pouched and are supplied multiple pouches per box.

IMPORTANT:

There are various techniques used in placing and affixing carotid bypass shunts. These shunts should be used only by qualified physicians thoroughly familiar with endarterectomy and patch-graft angioplasty techniques. Carotid bypass shunts are influenced by three significant factors:

1. **The size, configuration and disease state of the artery.**
2. **The design and dimensions of the shunt.**
3. **The stabilization technique to hold the shunt in place.**

Stabilization can be accomplished by many techniques, including vessel loops, umbilical-type tapes, heavy ligatures, tourniquets, or surgical clamps. Surgical clamps should be appropriate in size and design for the shunt chosen.

Warning:

When there is an incompatibility between the stabilizing technique, the artery and the shunt, slippage can occur. In order to avoid slippage and/or vessel damage or disruption when the clamp technique is utilized, users should ensure that the clamps they employ are compatible in size and design with the particular shunt they intend to use.

Suggested Insertion Techniques:**Non T-shaped Shunts**

The carotid arteries should be isolated, tourniquet loops placed, and the arterial incision made. The proper size shunt must be chosen, determined at the time of surgery, based on the surgeon's experience and training. Forcing a shunt that is too large into an artery may cause vessel disruption. Insert the shunt into the artery. Care should be taken not to dislodge plaque, clot, or to raise intimal flaps.

The bypass shunt is held in place by tension on the tourniquet loops. The endarterectomy is then performed. When the endarterectomy is completed, using standard cardiovascular techniques, which may include patch-graft angioplasty, the shunt is removed and the incision is closed.

T-shaped Shunts

The carotid arteries should be isolated, tourniquet loops placed, and the vessels occluded. A clamp is placed on the vertical limb of the T-shaped shunt. The artery is opened longitudinally. The proper size shunt must be chosen, determined at the time of surgery, based on the surgeon's experience and training. Forcing a shunt that is too large into an artery may cause vessel disruption. The distal end of the shunt is inserted into the internal carotid artery with a clamp placed on the proximal portion of the shunt. Allow blood from the distal internal carotid artery to fill the distal shunt and T-extension by briefly opening the clamps on the T-extension and the internal carotid artery.

Place a clamp on the distal portion of the shunt. Insert the proximal end of the shunt into the common carotid artery. Remove the clamp on the common carotid artery and the T-extension. Allow blood to fill the proximal portion of the shunt and T-extension. Confirm that the shunt is patent. Remove any air bubbles.

Replace the clamp on the T-extension. Remove the clamps on the distal shunt to allow the blood to flow through the shunt.

Warning: Assure that the shunt is properly stabilized in the artery.

Following the endarterectomy, the arteriotomy is closed, following standard cardiovascular techniques, and the shunt is removed.

Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid.

Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

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Bard, Javid, Brener, and Burbank are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.
Copyright ©2007 C. R. Bard, Inc. All Rights Reserved. Printed In The USA.
Copyright © 2007 C. R. Bard, Inc. Tous droits réservés. Imprimé aux Etats-Unis.
Copyright © 2007 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedrukt in den USA.
Copyright © 2007 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.
Copyright © 2007, C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.