<b>Instructions For Use</b>	
<b>Bard® Cardiovascular Fabrics</b> . . . . .	<b>1</b>
Mode d'emploi	
Tissus cardio-vasculaires Bard® . . . . .	3
Gebrauchsanweisung	
Bard® Gewebe für kardiovaskuläre Anwendungen . . . . .	5
Istruzioni per l'Uso	
Tessuti Cardiovascolari Bard® . . . . .	7
Instrucciones de uso	
Parches cardiovasculares Bard® . . . . .	9
Gebruiksaanwijzing	
Bard® cardiovasculaire weefsels . . . . .	11
Instruções de utilização	
Tecidos Cardiovasculares Bard® . . . . .	13
Οδηγίες χρήσης	
Καρδιαγγειακά Υφάσματα Bard® . . . . .	15
Brugervejledning	
Bard® Kardiovaskulære Stofstykker . . . . .	17
Bruksanvisning	
Bard® Kardiovaskulära textilier . . . . .	19
Käyttöohjeet	
Kardiovaskulaariset Bard®-kankaat . . . . .	21
Bruksanvisning	
Bard® kardiovaskulære stoffer . . . . .	23
Instrukcja użycia	
Tkaniny Bard® do użycia w chirurgii sercowo-naczyniowej . . . . .	25
Használati utasítás	
Bard® kardiovaskuláris foltok . . . . .	27
Pokyny k použití	
Kardiovaskulární látky Bard® . . . . .	29
Kullanma Talimatı	
Bard® Kardiyovasküler Kumaşlar . . . . .	31

## Bard® Cardiovascular Fabrics

### Instructions For Use

**Caution: Federal (U.S.A.) law restricts these devices to sale by or on the order of a physician. Only physicians qualified in cardiovascular techniques should use these products.**

#### Product Description:

**Bard** Cardiovascular Fabrics are available in various sizes, materials and types of construction.

The following Fabrics are constructed of **knitted polyester** (polyethylene terephthalate). Knitted construction is designed to resist fraying at the cut edges.

- **Bard *Sauvage***® Filamentous Fabric

Both sides of this fabric are filamentous, one considerably more than the other. The more defined filamentous surface is indicated by reference markings. These markings are at 2cm increments and may be used as a guide in cutting the fabric (see the Insertion Techniques).

- **Bard *DeBakey***® Double Velour Fabric

Both sides of this fabric contain a velour surface which is designed to enhance tissue ingrowth.

- **Bard *DeBakey***® Elastic Knit Fabric

This fabric is designed to stretch primarily in one direction and has high permeability.

The following fabric is constructed of woven polyester.

- **Bard *DeBakey***® Woven Fabric

This low permeability fabric has orientation markings on one side to assist in cutting the fabric to the desired size and shape (see Warning #5).

The following fabric is constructed of **knitted PTFE** (polytetrafluoroethylene).

- **Bard** Edwards Outflow Tract Fabric

This fabric is knitted with PTFE fibers and has an additional dispersion of PTFE to reduce the permeability of the fabric.

#### Indications for Use:

1. The knitted and woven polyester fabrics (**Bard *Sauvage*** Filamentous fabrics, **Bard *DeBakey*** Double Velour fabrics, **Bard *DeBakey*** Elastic Knit fabrics, **Bard *DeBakey*** Woven fabrics) are indicated for use in cardiovascular surgical procedures requiring patch graft angioplasty such as carotid endarterectomy. These fabrics are also indicated for repair of certain intracardiac anomalies such as septal defects.
2. **DeBakey** Woven Fabrics and Edwards Outflow Tract Fabrics are indicated for aortic and pulmonary outflow tract repair and repair of intracardiac defects.

#### Contraindications:

Due to the high permeability of the knitted polyester fabrics, they are contraindicated for use in patients requiring prolonged systemic or high dose heparinization except for use in repair of intracardiac defects.

#### Warnings:

1. These fabrics must be properly preclotted with non-heparinized blood prior to exposure to full arterial pressure to avoid unnecessary bleeding or blood loss. Preclotting is not required when the fabrics are used in the intracardiac position.
2. While it is common practice to employ moderate doses of intraoperative heparin, care should be taken not to exceed the manufacturer's recommended dose for such procedures. Excessive amounts of heparin can result in bleeding.
3. Due to variability in patient response to heparinization, it is essential that the adequacy of anticoagulation during surgery and the precision of neutralization at the conclusion of the procedure be closely monitored. Adherence to a strict protocol, determined by each hospital, can prevent excessive bleeding.
4. As with any cardiovascular fabric, occasional difficulties with hemostasis may occur. In the event that hemostasis is not easily obtained, the surgeon may wish to consider the following:
  - Investigate for systemic coagulopathy and treat appropriately.
  - Compression, as necessary.
  - Additional sutures and/or pledgets at the anastomosis, as necessary.
  - Reversal of heparin with protamine sulfate, as necessary.
  - Utilization of topical coagulation therapy, e.g., thrombin, as necessary.
5. As with all woven fabric constructions, care must be taken when trimming the fabric to minimize the potential for fraying at cut edges. Cautery is highly recommended for heat sealing on all woven patch edges. If the edges of the fabric are not heat sealed, then sutures must be at least 2mm from the cut edge.
6. Do NOT expose PTFE Fabrics to temperatures greater than 500°F (260°C). PTFE decomposes at elevated temperatures, producing highly toxic decomposition products.<sup>1</sup>

#### Precautions:

1. STERILE, unless the package is opened or damaged. Single use only.
2. These products are sterilized by ethylene oxide. DO NOT RESTERILIZE.
3. Care should be taken when clamping the fabrics to avoid damage to the fibers and preclot. Only atraumatic, shod-type vascular clamps should be used.
4. Taper point, noncutting needles are recommended so the fabric fibers will not be cut during suturing.
5. After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

#### Adverse Reactions:

Adverse reactions that may occur with the use of these products or with any cardiovascular implant procedure include perioperative hemorrhage, implant bleeding, tissue erosion, anastomotic aneurysms, and infection.

#### Equipment Required:

Scissors, suture, forceps, cautery, suture with taper point noncutting needles, and atraumatic shod-type vascular clamps.

## Directions For Use

#### How Supplied/Stored:

**Bard** Cardiovascular Fabrics are supplied sterile, unless package is damaged or opened. Single use, do not resterilize. **Bard** Cardiovascular Fabrics are single-pouched and are supplied multiple pouches per box. Refer to Warning #6 for storage conditions.

#### Preclotting Instructions:

1. Proper preclotting with non-heparinized blood is necessary to prevent bleeding through the fabric wall.
2. If the patient has been heparinized prior to preclotting the fabric, draw blood for the preclot and counteract the heparin by adding sufficient topical

- thrombin to initiate clotting. The fabric will then preclot readily.
- If bleeding is noted after the fabric has been preclotted, see Warnings for additional hemostasis information.
  - Suggested methods of preclotting are described below, however, it is appropriate to employ whatever your institution has accepted and validated as standard practice for permeable knitted and woven cardiovascular "fabrics".

#### **Sauvage® Filamentous Fabric**

Preclotting should be accomplished in the following 2-step manner:

- Place the patch in a basin and cover it with non-heparinized blood. Leave the patch submerged until the blood has clotted solidly. If, for some reason, the blood does not clot after six to seven minutes, add sufficient topical thrombin to cause solid clotting. Then take the patch from the pan and remove the visible clot from the surface.
- Place this patch (which has now been converted from a porous fabric to a highly thrombogenic, impervious fabric coated with fibrin) in a second basin and cover it with high dose heparinized blood (400 units/cc). The heparin will combine with the thrombin and render the thrombin inert. The patch should be left submerged in the heparinized blood until ready for use.

#### **Other Fabrics**

Immerse the fabric in a small quantity of the patient's own non-heparinized blood. Wait until the blood has solidly clotted. Gently remove visible clot from the surface.

#### **Insertion Techniques:**

##### **Sauvage Filamentous Fabrics**

##### **Reference Markings:**

- The more defined filamentous surface of the fabric is indicated by reference markings. This surface should be placed outward in patch graft procedures and toward the right side of the heart when used as an intracardiac patch.
- The arrow on the reference markings indicates the direction of the natural curvature and greatest stretch of the fabric. The arrows should be oriented at right angles to the long axis of the artery during patch graft procedures.

##### **Patch Graft Procedures:**

- Systemic heparinization of the patient is recommended for patch graft angioplasty procedures.
- Trim the fabric with scissors to the desired dimensions for the angioplasty procedure.
- While the suture technique employed will depend upon the surgeon's preference, a few general guidelines are appropriate:
  - Special attention should be given to avoid narrowing the artery at either end of the patch graft. In small vessels, interrupted sutures may be advisable.
  - Each suture should be made sufficiently deep into the heart or vascular wall to avoid its tearing out. Each suture should be placed at a 90 degree angle to the wall, passing the needle from the intimal surface outwards and including the adventitial layer for strength.
  - The cardiovascular wall/patch graft junction should be smooth, with particular care to avoid turning the free edge of the prosthesis inward into the lumen.
- A 15 to 20 minute delay in neutralizing systemic heparin by protamine administration is recommended.

#### **DeBaKey® Knitted Fabrics**

##### **Patch Graft Procedures:**

- Trim the fabric with scissors to the desired dimensions for the angioplasty procedure.
- The anastomosis of the patch graft to the arterial wall is preferably performed by a continuous through-and-through suture of nonabsorbable suture material with a double-ended needle, by attaching the patch graft to one end of the opening in the arterial wall.
- One end of the suture is then used to anastomose the patch graft to one edge of the arterial wall for about half of its extent. The other end is then used to anastomose the patch graft to the other edge of the arterial wall and the remaining edges of the first side.
- Just before completion of the procedure and before tying the two ends of the suture, the small remaining opening between the patch graft and the arterial wall is held closed with tissue forceps and the proximal occluding clamp on the artery is released to permit one or two pulses of blood to flow into the patch graft. This procedure is continued intermittently about every 30 to 60 seconds until the interstices of the patch graft are completely sealed and no further bleeding takes place through the patch graft or suture line, usually within several minutes.

**NOTE:** If bleeding continues after several minutes, see Warnings for additional hemostasis information.

#### **DeBaKey Woven Fabrics**

- See Warning #5.
- Suturing Techniques: See Insertion Techniques, **Sauvage** Filamentous Fabric, Patch Graft Procedures #3.

#### **WARRANTY**

Bard Peripheral Vascular warrants to the first purchaser of this product that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some states/countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your state/country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet. In the event that 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

#### **Reference:**

- Guide to the Safe Handling of Fluoropolymer Resins, 3rd Edition, The Fluoropolymers Division of the Society of the Plastics Industry, Inc.

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**Kaynakça:**

1. Guide to the Safe Handling of Fluoropolymer Resins, 3rd Edition, The Fluoropolymers Division of the Society of the Plastics Industry, Inc.

**Bard Sauvage Filamentous Knitted Polyester Fabric**

Tissu filamenteux Bard Sauvage en polyester tricoté  
Bard Sauvage  
Gestricktes Polyesterfasermaterial  
Tessuto filamentoso in maglia di poliester  
Bard Sauvage  
Parche de punto de poliéster filamentoso Bard Sauvage  
Bard Sauvage filamenteus gebreid polyesterweefsel  
Tecido filamentoso de malha de poliéster Bard Sauvage  
Νηματούδες πλεκτό ύφασμα Sauvage της Bard από πολυεστέρα  
Bard Sauvage strikket stofstykke af polyester med trådmønster  
Bard Sauvage fintrådig stickad polyestertextil  
Kuitumainen Bard Sauvage -polyesterineuleakangas  
Bard Sauvage trådmønstrede strikket polyesterstoff  
Włókienkowa dzianina poliesterowa Sauvage firmy Bard  
Bard Sauvage szálas kötött poliészter folt  
Vláknitá pletená polyesterová látka Bard Sauvage  
Bard Sauvage |pliksi  
Örme Polyester Kumaş

**Bard DeBaKey Elastic Knitted Polyester Fabric**

Tissu élastique Bard DeBaKey en polyester tricoté  
Bard DeBaKey  
Elastisches gestricktes Polyestermaterial  
Tessuto in maglia di poliester elastizzata  
Bard DeBaKey  
Parche de punto de poliéster elástico Bard DeBaKey  
Bard DeBaKey elastisch gebreid polyesterweefsel  
Tecido de malha elástica de poliéster Bard DeBaKey  
Ελαστικό πλεκτό ύφασμα DeBaKey της Bard από πολυεστέρα  
Bard DeBaKey elastisk strikket stofstykke af polyester  
Bard DeBaKey elastisk stickad polyestertextil  
Bard DeBaKey -polyesterjoustoneuleakangas  
Bard DeBaKey elastisk strikket polyesterstoff  
Elastyczna dzianina poliesterowa DeBaKey firmy Bard  
Bard DeBaKey rugalmas kötött poliészter folt  
Elastická pletená polyesterová látka Bard DeBaKey  
Bard DeBaKey Elastik  
Örme Polyester Kumaş

**Bard Edwards Outflow Tract Knitted PTFE Fabric**

Tissu pour voie d'éjection Bard Edwards en PTFE tricoté  
Bard Edwards  
Gestricktes PTFE-Material für Abflussgänge  
Tessuto in maglia di PTFE per tratto di efferenza Bard Edwards  
Parche de punto de PTFE para vías de salida Bard Edwards  
Bard Edwards  
gebreid PTFE weefsel voor outflowtraject  
Tecido das vias de ejección de malha de PTFE Bard Edwards  
Πλεκτό ύφασμα χώρου εξώθησης Edwards της Bard από PTFE  
Bard Edwards strikket PTFE-stofstykke til udløbssystem  
Bard Edwards textil för utflödesbanor av stickad PTFE  
Bard Edwards Outflow Tract -PTFE-neuleakangas  
Bard Edwards Outflow Tract Strikket PTFE-stoff  
Dzianina PTFE do dróg odpływowych Edwards firmy Bard  
Bard Edwards kötött PTFE kiáramlási pálya folt  
Pletená PTFE látka Bard Edwards pro oblast výtoku  
Bard Edwards Çıkış Kanalı  
Örme PTFE Kumaş

**Bard DeBaKey Double Velour Knitted Polyester Fabric**

Tissu de velours double Bard DeBaKey en polyester tricoté  
Bard DeBaKey  
Gestricktes Doppelveloursmaterial aus Polyester  
Tessuto in maglia di poliester doppio velour  
Bard DeBaKey  
Parche de punto de poliéster de doble terciopelo Bard DeBaKey  
Bard DeBaKey dubbel velours gebreid polyesterweefsel  
Tecido de veludo duplo de malha de poliéster Bard DeBaKey  
Πλεκτό ύφασμα από διπλό βελούρο DeBaKey της Bard από πολυεστέρα  
Bard DeBaKey dobbelt velour strikket stofstykke af polyester  
Bard DeBaKey dubbel velourstickad polyestertextil  
Bard DeBaKey -kaksoisveluuri-polyesterineuleakangas  
Bard DeBaKey dobbelt velurstrikket polyesterstoff  
Podwójna welurowa dzianina poliesterowa DeBaKey firmy Bard  
Bard DeBaKey kétrétegű velúr kötött poliészter folt  
Dvojitá velurová pletená polyesterová látka Bard DeBaKey  
Bard DeBaKey Çift Katlı Velur  
Örme Polyester Kumaş

**Bard DeBaKey Woven Polyester Fabric**

Tissu Bard DeBaKey en polyester tissé  
Bard DeBaKey  
Polyestergewebe  
Tessuto in poliester Bard DeBaKey  
Parche de poliéster entretejido Bard DeBaKey  
Bard DeBaKey  
geweven polyestertextiel  
Tecido de poliéster tecido Bard DeBaKey  
Υφαντό ύφασμα DeBaKey της Bard από πολυεστέρα  
Bard DeBaKey  
vævet stofstykke af polyester  
Bard DeBaKey  
vävd polyestertextil  
Bard DeBaKey -polyesterikudoskangas  
Bard DeBaKey vevd polyesterstoff  
Tkanina poliesterowa DeBaKey firmy Bard  
Bard DeBaKey  
szőtt poliészter folt  
Tkaná polyesterová látka Bard DeBaKey  
Bard DeBaKey  
Dokuma Polyester Kumaş

**Attention, See Instructions For Use**

Attention, voir le mode d'emploi  
Achtung! Siehe Gebrauchsanweisung  
Attenzione, leggere le istruzioni per l'uso  
Atención: consultar las instrucciones de uso  
Let op, lees de gebruiksaanwijzing  
Atenção, consulte as Instruções de Utilização  
Προσοχή, βλ. Οδηγίες χρήσης  
Bemærk, Se brugervejledning  
Obs! Se bruksanvisningen  
Huomio! Lue käyttöohjeet  
OBS! Se bruksanvisningen  
Uwaga: Należy zapoznać się z instrukcją użytkowania  
Figyelem, lásd a használati utasítást  
Pozor, viz pokyny k použití  
Dikkat, Kullanma Talimatına Bakınız

**Lot Number**

Numéro du lot  
Chargennr.  
Numero di lotto  
Número de lote  
Serienummer  
Número de Lote  
Αριθμός Παρτίδας  
Lot-number  
Produktionsnr  
Eränumero  
Lot nummer  
Numer serii  
Tételszám  
Číslo šarže  
Parti Numarası

**Single Use**

A usage unique  
Nur zur einmaligen Verwendung  
Monouso  
Para un solo uso  
Eenmalig gebruik  
Utilização única  
Για μία μόνο χρήση  
Engangsbrug  
För engångsbruk  
Vain kertakäyttöön  
Til engangs bruk  
Do jednorazowego użytku  
Egyszeri használatra  
Jednorázové použití  
Tek Kullanımiktir

**Contents**

Contenu  
Inhalt  
Contenuto  
Contenido  
Inhoud  
Conteúdo  
Περιεχόμενα  
Inndhold  
Innehåll  
Sisällys  
Innhold  
Zawartość  
Tartalom  
Obsah  
İçerikler

**Do Not Use If Package Is Damaged Or Opened**

Ne pas utiliser si l'emballage est ouvert ou endommagé  
Nicht verwenden, wenn die Verpackung geöffnet oder beschädigt ist.  
Non utilizzare se la confezione è aperta o danneggiata  
No usar si el envase está abierto o dañado.  
Niet gebruiken wanneer de verpakking geopend of beschadigd is.  
Nào utilizez se a embalagem estiver aberta ou danificada.

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**Copyright ©2008 C. R. Bard, Inc. All Rights Reserved. Printed In The USA.**

Copyright 2008 C. R. Bard, Inc. Tous droits réservés. Imprimé aux Etats-Unis.  
Copyright ©2008 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.  
Copyright ©2008 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.  
Copyright ©2008 C. R. Bard, Inc. Reservados todos los derechos. Impreso en EE.UU.  
Copyright ©2008 C. R. Bard, Inc. Alle rechten voorbehouden. Gedrukt in de VS.  
Copyright ©2008 C. R. Bard, Inc. Todos os direitos reservados. Impresso nos EUA.  
Πνευματικά δικαιώματα ©2008 C. R. Bard, Inc. Με την επιφύλαξη παντός δικαιώματος. Εκτυπώθηκε στις ΗΠΑ.  
Copyright ©2008 C. R. Bard, Inc. Alle rettigheder forbeholdes. Trykt i USA.  
Copyright ©2008 C. R. Bard, Inc. Eftertryck förbjudes. Tryckt i USA.  
Copyright ©2008 C. R. Bard, Inc. Kaikki oikeudet pidätetään. Painettu USA:ssa.  
Copyright ©2008 C. R. Bard, Inc. Alle rettigheter förbehållt. Trykt i USA.  
Copyright ©2008 C. R. Bard, Inc. Wszelkie prawa zastrzeżone. Wydrukowano w USA.  
Szerzői jog ©2008 C. R. Bard, Inc. Minden jog fenntartva. Nyomtatva az USA-ban.  
Copyright ©2008 C. R. Bard, Inc. Veškerá práva vyhrazena. Vytiskněno v USA.  
Telif Hakkı ©2008 C. R. Bard, Inc. Tüm Hakları Saklıdır. A.B.D.'de Basılmıştır.

**Sterilized By Using**

**Ethylene Oxide**  
Stérilisé à l'oxyde d'éthylène  
Sterilisiert mit Ethylenoxid  
Sterilizzato mediante ossido di etilene  
Sterilizado usando óxido de etileno  
Gesteriliseerd met behulp van ethyleenoxide  
Esterilizado por óxido de etileno  
Αποστειρωμένο με αιθυλενοξείδιο  
Steriliseret ved ætylenoxid  
Steriliserad med etylenoxid  
Steriloitu etyleenioksidilla  
Sterilisiert med etylenoksid  
Produkt sterylizowany tlenkiem etylenu  
Etilénoxidál sterilizálva  
Sterilizovano etylenoxidem  
Etilen Oksit ile Sterilize Edilmiştir

**E.E.A. Authorized Representative**

Représentant autorisé E.E.E.  
EEA Berechtigter Vertreter  
Rappresentante autorizzato EEA  
Representante autorizado en el EEE  
Gemachtigde binnen de EER  
Representante autorizado no EEE  
Εξουσιοδοτημένος αντιπρόσωπος στον ΕΟΧ  
EÖF-autoriseret repræsentant  
Auktoriserad representant inom EES  
ETA:n valtuutettu edustaja  
EØS-autoriseret representant  
Autoryzowany przedstawiciel dla krajów Europejskiego Obszaru Gospodarczego  
EGT meghatalmazott képviselő  
Autorizovaný zástupce pro EHP  
Avrupa Ekonomik Alanı Yetkili Temsilcisi

**Bard and DeBakey are registered trademarks of C. R. Bard, Inc. or an affiliate. Sauvage is a registered trademark licensed to C. R. Bard, Inc. or an affiliate.**

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**Bard Peripheral Vascular, Inc.**
Bard Limited