# CENTURION®

## PTA Balloon Dilatation Catheter

### ENGLISH

#### INSTRUCTIONS FOR USE

**CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.**

#### Device Description

The **Centurion** PTA Balloon Dilatation Catheter is a high pressure balloon catheter consisting of an over the wire catheter with a balloon fixed at the distal tip. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The catheter includes an atraumatic tip to facilitate advancement of the catheter to and through the stenosis. The over the wire, dual lumen catheter is compatible with 0.035" guidewire and is available in 50, 75, and 120 cm working lengths. The proximal portion of the catheter includes a female luer lock hub connected to the inflation lumen, and a female luer-lock hub connected to the guidewire lumen. Packaged with every product is a stylet and a profile reducing sheath that is positioned over the balloon for protection before use. This product is not manufactured with any latex.

#### Indications for Use

**Centurion** PTA Balloon Dilatation Catheters are recommended for Percutaneous Transluminal Angioplasty (PTA) of the renal, iliac, and femoral arteries, post-deployment dilatation of peripheral vascular stents<sup>1</sup>, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This catheter is not for use in coronary arteries.

#### Contraindications

None known.

#### Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). **Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess or re-sterilize.**
2. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
3. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
4. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over pressurization, use of a pressure monitoring device is recommended.
5. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

#### Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The **Centurion** PTA Balloon Dilatation Catheter shall only be used by physicians trained in the performance of percutaneous transluminal angioplasty.
3. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size sheath introducer than indicated on the label.
4. Use the recommended balloon inflation medium. Never use air or other gaseous medium to inflate the balloon.

5. If resistance is felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
6. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
7. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze and rinsed with sterile normal saline.

#### Potential Adverse Reactions

The complications which may result from a peripheral balloon dilatation procedure include:

- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Embolization
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

#### Directions for Use

##### Handling and Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources. Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date. Do not use if packaging is damaged or opened.

#### Equipment Required

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- 0.035" guidewire

#### Dilatation Catheter Preparation

**NOTE: Remove the stylet prior to removing the protective sheath as part of the preparation for use.**

1. Remove Catheter from package. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
2. Remove the balloon guard by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon protector and slide distally off of the balloon catheter.
3. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium. Do not use air or any gaseous medium to inflate the balloon.
4. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
5. Connect the syringe to the stopcock.
6. Hold the syringe with the nozzle pointing downward, open the stopcock and aspirate for approximately 15 seconds. Release the plunger.
7. Repeat step #6 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
8. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

#### Use of the **Centurion** PTA Balloon Dilatation Catheter

1. Backload the distal tip of the **Centurion** PTA Balloon Dilatation Catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over the wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized

dilatation catheter.

3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and inflate the balloon to the appropriate pressure.
4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
5. While maintaining negative pressure and the position of the guidewire, withdraw the deflated dilatation catheter over the wire through the introducer sheath. Use of a gentle counterclockwise motion may be used to help facilitate catheter removal through the introducer sheath.

#### Balloon Reinsertion

**Precaution:** Do not continue to use the balloon catheter if the shaft has been bent or kinked.

**Precaution:** Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze and rinsed with sterile normal saline.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
3. Continue the procedure according to the "Use of the **Centurion** PTA Balloon Dilatation Catheter" section herein.

**Warning: After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.**

#### Warranty

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product is not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

<sup>1</sup> All stents should be deployed in accordance with the manufacturer's indications and instruction for use. Laboratory test have demonstrated that Centurion PTA Balloon Dilatation Catheters have been shown to withstand multiple inflations to rated burst pressures within previously deployed Palmaz® Stents. Data on file.

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**Copyright © 2007 C. R. Bard, Inc. All rights reserved. Printed in the USA.**
Copyright © 2007 C. R. Bard, Inc. Alle rettigheter forbeholdt. Trykkt i USA.
**Bard and Centurion are registered trademarks of C.R. Bard, Inc. or an affiliate.**
**Palmaz is a registered trademark of Johnson & Johnson**
**Bard Peripheral Vascular, Inc.**
Bard Limited