# ***Vectra***<sup>®</sup> **Vascular Access Graft**

**English  
Français  
Deutsch  
Italiano  
Español**

## ***INSTRUCTIONS FOR USE***

14160.Q 11/02

## INSTRUCTIONS FOR USE


### Sterile (EO) – For Single Use Only

*Caution: Federal (USA) law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).*

### Description


The Vectra Vascular Access Graft (VAG) is composed of Thoralon®, a proprietary blend of segmented polyetherurethaneurea and a siloxane containing surface-modifying additive (SMA). The grafts are available in two diameters, 5.0mm or 6.0mm, each with a wall thickness of approximately 1.0mm.

The Vectra VAG features a self-bonded tri-layer design of Thoralon. (Figure 1)

**Figure 1. Cutaway view of Vectra VAG**

- The microporous, blood-contacting inner layer contains an SMA rich environment designed to minimize platelet adhesion.
- The solid middle layer prevents any communicating pathway between the inner and outer layers, and gives the graft its strength, flexibility, and self-sealing properties.
- The microporous outer layer is designed to promote tissue ingrowth for enhanced graft anchoring, and is reinforced for kink resistance with two thin monofilament polyester fibers. Each reinforcement monofilament is wound circumferentially in graduated concentrations along the length of the graft, with an unreinforced area at either end. The concentration of the wraps increases from 10-12 wraps per inch (W.P.I.) or 4-5 wraps per centimeter (W.P.C.) at either end of the graft to 24-32 W.P.I. or 9.5-12.5 W.P.C. in the mid-portion of the graft for added kink resistance in the apex when used in a loop configuration. (Figure 2)

**Figure 2. Vectra graft design (40cm pictured)**

### Indications for Use

The Vectra VAG is indicated for use as a subcutaneous arteriovenous conduit for blood access.

### Contraindications

None.

## Potential Complications with the Use of a Vascular Prosthesis

Potential complications that may occur with any surgical procedure involving a vascular prosthesis include, but are not limited to: aneurysm; anastomotic disruption or tearing of the suture line and/or host vessel; embolic events; infection; bleeding; occlusion; stenosis; thrombosis; kinking/compression; swelling of the implanted limb; formation of hematomas or pseudoaneurysms; steal syndrome; and/or skin erosion.

### Warnings

**AVOID EXCESSIVE AXIAL ELONGATION OR STRETCHING OF THE GRAFT (>10%) DURING HANDLING AT IMPLANTATION.** Cut the graft long enough to prevent stress on the anastomoses and/or allow for a full range of body movements. Excessive elongation or stretching of the graft may result in damage to the graft microporous layers or anastomotic disruption that could lead to hematoma, bleeding, pseudoaneurysm, ischemia, or loss of limb or limb function.

Do NOT use product if package has been damaged or opened, as sterility may be compromised.

Do NOT resterilize the graft.

### Precautions

1. When using a balloon angioplasty or embolectomy catheter within the lumen of the Vectra VAG, match the inflated balloon size to the inner diameter size of the VAG. Over-inflation of the balloon or use of an inappropriately sized balloon catheter may damage the graft. Care must be exercised to avoid causing excessive axial elongation of the graft during retraction.
2. Preclotting this graft is not necessary.
3. It is necessary to hydrate the graft in a sterile solution of normal physiologic saline prior to implantation. The inner and outer luminal walls of the graft are microporous and the voids in these surfaces contain air that must be displaced. See "Operative Techniques" for more information.
4. Use of the Vectra sheath tunneler is necessary to minimize subcutaneous trauma and to minimize the force required to position the graft during implant. Tunneling permits graft placement without pulling that can produce excessive force that could result in damage to the graft microporous layers.
5. If clamping is necessary, use only atraumatic or appropriate vascular smooth-jawed or shod clamps to avoid damage to the graft wall during implantation.
6. Because of the monofilament reinforcement, only an end-to-end, not an end-to-side, anastomosis of the revision graft to the implanted graft is recommended.
7. Guideline 17 from the Clinical Practice Guidelines from the Final Report from the Vascular Access Work Group of the National Kidney Foundation (US) – Dialysis Outcomes Quality Initiative, recommends the following for infected AV grafts: *"When infected, a dialysis graft should be treated surgically. An untreated access infection may lead to bacteremia, sepsis, hemorrhage, and death. Surgical exploration and removal of any infected graft or graft segment is necessary for resolution of the infection because the graft material acts as a foreign body unless eliminated."*

Patients receiving grafts in the femoral location are known to have more infections and thrombosis complications, and should be monitored for signs of infection and attended to promptly if any are identified.<sup>1</sup>

## Summary of the Vectra® VAG Clinical Experience

**Purpose:** The Thoratec Vectra VAG was compared to ePTFE grafts in a prospective randomized study to establish equivalency in the safety and effectiveness for the Vectra and ePTFE grafts when used as subcutaneous arteriovenous conduits.

---

1. Taylor et al. Results and complications of arteriovenous access dialysis grafts in the lower extremity: A five year review. Am. Surg. 1996 Mar;62(3):188-91.

**Study Design:** In a prospective, randomized controlled study, a total of 142 patients (*Vectra*-70, ePTFE-72) were enrolled at five institutions in the US. Patient characteristics between the two groups were matched with respect to risk factors and demographics. All grafts implanted in this study were 6mm in diameter.

Primary endpoints included primary patency and functional patencies (secondary patency), time to first thrombosis, and adverse events. In addition, the study also documented the time to first dialysis access and time to hemostasis after cannulation. Time to first thrombosis, and primary and secondary patencies were determined from the date of implant as well as the date of first dialysis access, because a significant number of *Vectra* grafts were cannulated earlier than the control group.

**Table 1. Primary & Secondary Patencies of *Vectra* and ePTFE grafts**

Characteristic	<i>Vectra</i> (365 days)	ePTFE (365 days)	p-value
Survival to first thrombosis			
- from implant date	42%	50%	0.14
- from first dialysis date	43%	47%	0.46
Primary patency			
- from implant date	36%	36%	0.33
- from first dialysis date	39%	30%	0.77
Functional patency			
- from implant date	75%	80%	0.55
- from first dialysis date	83%	82%	0.75

As shown in Table 1, the primary and functional patencies for *Vectra* grafts are substantially equivalent to those for ePTFE grafts for both implant date and first dialysis date outcomes. Major complications for both *Vectra* and ePTFE grafts are presented in Table 2.

**Table 2. Rate of Major Complications During the Study**

Complication	<i>Vectra</i>		ePTFE	
	# events	(#Pts)	# events	(#Pts)
<b>General Events</b>				
Anastomotic Obstruction	51	(27)	35	(23)
Aneurysm	0	(0)	0	(0)
Infection, Surgical	5	(4)	6	(6)
Intra-graft Obstruction	13	(11)	17	(15)
Kinking	7	(7)	0	(0)
Pseudoaneurysm	0	(0)	3	(3)
Wound Healing/Skin Erosion	2	(2)	0	(0)
Sub Total Complications <sup>1</sup>	78	(37)	61	(35)
<b>Native Vessel Disease</b>				
Arterial Disease	6	(6)	2	(2)
<b>Central and Distal Vein Disease</b>				
Non-vascular arm related <sup>2</sup>	10	(8)	3	(2)
<b>Virginal arm incidents</b>				
Pre-existing Conditions <sup>3</sup>	28 <sup>4</sup>	(13)	12	(8)
No Reported Conditions	10	(7)	13	(10)
Sub Total Complications	54	(34)	30	(22)
<b>Patient Events<sup>5</sup></b>	7	(6)	2	(2)
<b>Deaths</b>	11	(11)	10	(10)
<b>Technical events<sup>6</sup></b>	18	(15)	8	(7)

1. Patient within category may have had more than one complication.
2. Non-vascular arm is defined as prior implantation of a PTFE graft or the multiple placements of fistulas (indicating multiple failed fistulas) in the same arm as the study graft.
3. Stenotic disease indicated by prior history or diagnostic procedures, temporary venous catheter placement in the same side as study graft or history of multiple TVC placements, and history of multiple prior shunt failures.

4. Five complicated patients experienced 16 events due to multiple underlying patient factors.
5. Includes: edema, hypotension, hypercoagulable condition, and ecchymosis.
6. Includes compression (external), hematoma, kinking due to technique, steal syndrome, laceration, and nerve release procedure.

**Table 3. Early Access for Cannulation**

Characteristic	<b>Vectra N=63/70</b>	<b>ePTFE N=61/72</b>	<b>p-value</b>
Grafts accessed within			
- 3 days	33.3%	0.0%	<0.05
- 8 days	54.0%	0.0%	<0.05
- 14 days	68.3%	9.8%	<0.05

**Table 4. Time to Hemostasis after Cannulation**

Characteristic	<b>Vectra N=1475</b>	<b>ePTFE N=1341</b>	<b>p-value</b>
Percent of grafts with time to hemostasis			
- ≤ 2 minutes	Arterial – 44.7% Venous – 52.2%	Arterial – 11.8% Venous – 13.1%	<0.05 <0.05
- ≤ 5 minutes	Arterial – 83.5% Venous – 86.5%	Arterial – 34.7% Venous – 39.5%	<0.05 <0.05

There were a significantly higher number of *Vectra* group patients who had early use of the graft, than in the ePTFE group (Table 3). After cannulation, there were more patients in the *Vectra* group that had hemostasis times ≤ 2 minutes (Table 4). This benefit was enhanced when the time to hemostasis was ≤ 5 minutes. In spite of this early access, primary and functional patencies were similar in both *Vectra* and ePTFE groups.<sup>1</sup>

## Operative Techniques

### Opening the Package

Holding the outer tray, peel back the lid and remove the inner tray. Peel back the inner tray lid and carefully remove the graft using sterile atraumatic instruments or gloves.

### Hydrating the Graft

Prior to implantation, hydrate or soak the graft in a sterile solution of normal physiologic saline. (Figure 3) While it is soaking, gently compress the submerged graft to displace the air from the voids in the microporous structure. When the bubbles stop, the air has been displaced and you may proceed with graft placement.

**Figure 3. Hydrating the Vectra graft**

1. Data on File.

### Implantation Tips and Techniques


Due to the compliant nature of the *Vectra*, it will expand approximately 0.5cm once exposed to normal arterial pressure. We recommend the following tips and techniques during the implantation procedure in order to allow for graft expansion and reduce the incidence of kinking. These are further detailed on the following pages 7 through 11.

- The graft features a longitudinal orientation line, and a dashed line and apex marker to identify the more densely reinforced area, to facilitate proper graft positioning. (Figure 4)

**Figure 4. Apex marker**

- When performing the anastomoses, trim the graft to a length approximately 0.5 to 1.0cm shorter than the length measured. Pull each end slightly taut when performing anastomosis to vessel. The end of the graft should be beveled to accommodate a smooth lie of the graft.
- When positioning and trimming the graft, avoid bending those areas of less dense reinforcement, especially near the anastomotic ends. The non-reinforced segment must be trimmed to allow for proper sizing of the graft.
- Liberally irrigate with saline the inside of the tunneler sheath to facilitate freely pushing the graft through the sheath. (Figure 5)

**Figure 5. Lubricating tunneler sheath**

#### For looped configurations:

- Use of a sterile surgical skin marker to demarcate the proposed graft layout on the arm will facilitate final positioning of the graft.
- Use of the two-tunneler sheath technique facilitates movement of the graft for proper positioning prior to removing the sheath. (Figures 6 – 11)
- Perform one horizontal incision at the volar site, distal to the intended loop position, and dissect approximately 2cm of extra space in the subcutaneous pocket around the apex of the looped portion of the graft to allow for free expansion of the graft after establishing flow. (Figure 10)
- Use the apex marker to properly position the more densely reinforced portion (24-32 W.P.I./9.5-12.5 W.P.C.) of the graft at the loop apex. (Figure 9)
- We recommend a loop diameter of at least 3.5cm to minimize the incidence of kinking in the loop region.

#### Tunneling

Because the graft can be damaged if pulled excessively, it is necessary to use a *Vectra* VAG sheath tunneler with this graft. The tunneler is designed to permit graft placement by gently pushing the graft through the tunneler sheath rather than pulling. The *Vectra* VAG tunnelers (except for the full curve tunneler which can only be used with the 5mm) can be used with either the 5mm or 6mm *Vectra* VAG. Refer to the *Vectra* Tunneler IFU.

Liberally irrigate the inside of the tunneler sheath as well as all surfaces of the graft with a sterile solution of normal physiologic saline to facilitate slipping the graft through the tunneler (also see “Precautions” statement #3).

Always create tunnels at appropriate depths that will allow easy visualization and palpation of the Vectra VAG.

**Care should be taken to push rather than pull the graft through the tunnel.**

The Vectra VAG can be implanted using either a one-sheath or two-sheath technique as described below.

**One-sheath Technique:**

1. Place the sheath over the rod and thread the tunneler tip onto the end of the tunneler rod to hold the sheath in place.
2. Make two (2) incisions for distal and proximal entries of the implant site.
3. Insert the fully assembled tunneler into one of the incisions to create a tunnel subcutaneously between the distal and proximal incisions.
4. With the tip of the tunneler exposed, the tip is removed and the rod is pulled away, leaving the sheath in place subcutaneously.
5. Insert the Vectra VAG into the subcutaneously placed sheath and gently push the graft through using a rotating motion if necessary, while irrigating liberally with sterile saline.

**Do NOT pull the graft through the tunneler sheath. Do NOT twist the graft.**

6. For a looped graft configuration, position graft to make sure the loop portion is smooth and does not kink and that the more densely reinforced portion of the graft is at the loop apex. The center of the apex is marked as pictured in Figure 4.

**Do not attempt to reposition the graft after sheath removal.**

7. When positioning of graft is complete, carefully remove the sheath, leaving the graft in place subcutaneously.
8. Repeat steps 3 through 6 if using a looped graft configuration.

**Two-sheath Technique for Looped Graft Position:**


The following two-sheath technique has been successfully used with the Vectra VAG to facilitate the final positioning of the graft once the tunneler sheaths are in place in the subcutaneous tissue. (Figure 6)

**Figure 6. Two-sheath implantation technique**

This technique requires the use of a fully assembled tunneler and an extra sheath.

1. Place the sheath over the rod and thread the tunneler tip onto the end of the tunneler rod to hold the sheath in place.
2. Make two (2) incisions for distal and proximal entries of the implant site.
3. Insert the fully assembled tunneler into one of the incisions to create a tunnel subcutaneously between the distal and proximal incisions. (Figure 7)


**Figure 7. Insertion of fully assembled tunneler**

4. With the tip of the tunneler exposed, the tip is removed and the rod is pulled away, leaving the sheath in place subcutaneously.
5. Place the extra sheath over the rod and thread the tunneler tip onto the end of the tunneler rod to hold the sheath in place.
6. Repeat steps 3 and 4 for the remaining leg of the loop.
7. Insert the *Vectra* VAG into one side of the subcutaneously placed sheath and gently push the graft through while using a rotating motion if necessary, while irrigating liberally with sterile saline.


**Do NOT pull the graft through the tunneler sheath. Do NOT twist the graft.**

8. Leaving the first sheath cited in step 7 in place, insert the *Vectra* VAG into the other subcutaneously placed sheath and gently push the graft through while irrigating liberally. (Figure 8)

**Figure 8. Insert graft into lubricated sheath**

9. Position the graft to make sure the loop portion is smooth and does not kink and that the more densely reinforced portion (as shown in figure 4) of the graft is at the loop apex. (Figure 9)


**Figure 9. Positioning graft apex**

**Note:** Due to the compliant nature of the *Vectra*, it will expand slightly once exposed to normal arterial pressure. Therefore, to avoid kinking after pressurization, it is recommended that extra space of approximately 2cm be created in the subcutaneous pocket around the apex of the looped portion of the graft. Once the graft is pressurized, this will allow for free expansion of the graft.

- Retract the sheaths while gently pushing the graft loop into the pocket, again assuring the loop portion is smooth and that the more densely reinforced area is positioned at the loop apex. (Figure 10)


**Figure 10. Graft loop in pocket**

- Carefully remove the first sheath followed by the other sheath, leaving the graft in place subcutaneously. (Figure 11)

**Do not attempt to reposition the graft after sheath removal.**

**Figure 11. Sheath removal**


- After sheath removal, do not pull the graft or attempt to reposition.

**The Vectra Tunneler is a reusable device. Clean and sterilize before each use using the procedures listed in the Vectra Tunneler Instructions for Use.**

### **Anastomotic Preparation**

Once the graft is properly positioned, an anastomosis should be made between the vessel and the graft. The end of the graft should be beveled to accommodate a smooth lie of the graft. The anastomosis can be performed using a one or two-suture technique. The other end of the graft should be trimmed and anastomosed in the same manner. Note that the graft may be cut and/or sutured at any point along its reinforced length (Figure 12). (Also see "Precautions" statement #5.)

**Figure 12. Completed graft anastomosis**

### **Suturing**

Best results are achieved using a tapered, noncutting needle with a nonabsorbable monofilament 5.0 or 6.0 suture approximately the same size as the needle. Care must be taken to follow the curve of the needle and pull the suture at a 90° angle assuring the needle penetrates all three layers of the graft to minimize suture hole elongation and bleeding. Use systemic or local heparinization unless contraindicated.

Avoid excessive tension that may cause suture holes to elongate or tear; undue tension on the suture line; gaps between the graft and host vessel; and inappropriate

suture placement and bites. Elongation of suture holes or gaps between the graft and host vessel could lead to anastomotic bleeding.

## Thrombectomy

The *Vectra* VAG may be declotted as follows:

- Follow the catheter manufacturer's instructions regarding size, selection, and balloon inflation, matching the balloon size to the inner diameter of the graft. Over-inflation and excessive pulling may dilate or damage the graft.
- If devices such as Adherent Clot Catheter are used to declot the *Vectra* grafts, it is important to match the catheter size to the internal diameter of the graft.
- If a longitudinal incision is used, place stay stitches at each end of the incision before introducing the embolectomy catheter.
- If a transverse incision is used, no stay stitch is necessary and a horizontal mattress suture technique will aid in closure.

**Note:** Do not place undue stress on the anastomosis or incision when placing or removing the catheter.

## Surgical Revision

If it is necessary to repair the *Vectra* graft with a surgical interposition bypass graft, the use of a *Vectra* revision graft with uniform reinforcement is recommended. Two types of revision grafts are available in either 5mm or 6mm diameters:

- 10-12 W.P.I. (4-5 W.P.C.) reinforcement for the portion of the grafts with the less reinforced section, generally near the access sites.
- 24-32 W.P.I. (9.5-12.5 W.P.C.) for the portion of the graft with more densely reinforced section, generally in the loop portion of graft.

**Because of the monofilament reinforcement, only an end-to-end, not an end-to-side, anastomosis of the revision graft to the implanted graft is recommended.**

## Cannulation

**Note:** The *Vectra* Vascular Access Graft may be punctured for vascular access within 24 hours after implant, providing no contraindications are present. (see #4 below)

Insert the blood access (dialysis) needle at a 45° angle with the bevel up until the graft is penetrated. If the blood access needle is inserted such that the angle between the needle axis and the graft is too small, tears in the wall of the graft can occur. If the needle is inserted at a 90° angle, it increases the possibility of puncturing the far wall of the graft, which may lead to hematoma formation.

For best results follow the established cannulation practices listed below:

1. Rotate Cannulation Sites. Repeated cannulation in the same area may lead to damage of the graft wall and/or formation of hematoma or pseudoaneurysm. Needle puncture sites should be equally spaced along the subcutaneous length of the graft.
2. Do not cannulate within the dialysis needle's length of the proximal or distal anastomoses.
3. Strict adherence to aseptic technique is required to minimize infection.
4. As with all dialysis, do not cannulate if there are any signs of infection, bleeding, swelling, edema, hematoma, or in the absence of a strong "thrill".

After needle withdrawal, use **gentle, non-occlusive digital pressure** to compress the cannulation site to aid in hemostasis. The *Vectra* VAG seals quickly (1-5 minutes). Prolonged compression or use of stasis clamps may lead to clot formation, restricting flow through the graft.

Please refer to the *Vectra* 7 Vital Intervention Principles (For Dialysis) brochure #D028-1201 for additional information.

## **Storage**

To provide maximum protection, store the grafts in their original, unopened packages at room temperature. Avoid excessive heat or cold (>-10°C to <50°C).

The grafts must be used before the expiration date printed on each label.

## **How Supplied**

*Vectra* VAGs are available in either 5mm or 6mm diameters. All *Vectra* VAGs are supplied sterile.

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Thoratec, Thoralon, the Thoratec logo and Vectra are registered trademarks of Thoratec Corporation. Bard and IMPRA are registered trademarks of C. R. Bard, Inc., or an affiliate. Printed in USA. Part no. 14160.Q 11/02
Thoratec Corporation
IMPRA, Inc.
A Subsidiary of C. R. Bard, Inc.
Thoratec Europe Limited