



Conquest† PTA Catheter

Ordering Information for Canada

WESTERN CANADA	Tel. 800 268 2862	Fax 1-800-632-2109
ONTARIO	Tel. 800 387 9473	Fax 1-800-632-2109
EASTERN CANADA	Tel. 800 387 7851	Fax 1-800-632-2109

Shaft Length		Baloon Size		RBP** (atm)	Nominal Pressure* (atm)	Sheath Size (Fr)	Pack Quantity
50 cm	75 cm	Diameter (mm)	Length (cm)				
CQ-5052	CQ-7552	5	2	30	8	6	1
CQ-5054	CQ-7554	5	4	30	8	6	1
CQ-5058	CQ-7558	5	8	27	8	6	1
CQ-5062	CQ-7562	6	2	30	8	6	1
CQ-5064	CQ-7564	6	4	30	8	6	1
CQ-5068	CQ-7568	6	8	27	8	6	1
CQ-5072	CQ-7572	7	2	30	8	6	1
CQ-5074	CQ-7574	7	4	30	8	6	1
CQ-5076	CQ-7576	7	6	27	8	6	1
CQ-5078	CQ-7578	7	8	27	8	6	1
CQ-5082	CQ-7582	8	2	27	8	6	1
CQ-5083	CQ-7583	8	3	27	8	6	1
CQ-5084	CQ-7584	8	4	27	8	6	1
CQ-5086	CQ-7586	8	6	25	8	6	1
CQ-5088	CQ-7588	8	8	25	8	6	1
CQ-5092	CQ-7592	9	2	26	8	7	1
CQ-5094	CQ-7594	9	4	26	8	7	1
	CQ-75102	10	2	24	8	7	1
	CQ-75104	10	4	24	8	7	1
	CQ-75122	12	2	20	8	8	1
	CQ-75124	12	4	20	8	8	1

* Nominal Pressure: the pressure at which the balloon reaches its labelled diameter.

** RBP (Rated Burst Pressure): the pressure at which BARD has 95% confidence that 99.9% of the balloons will not burst at or below upon single inflation.

† Bard and Conquest are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.