



Bard† DeBakey† Double Velour Fabric

Ordering Information for Canada

WESTERN CANADA	Tel. 800 268 2862	Fax 1-800-632-2109
ONTARIO	Tel. 800 387 9473	Fax 1-800-632-2109
EASTERN CANADA	Tel. 800 387 7851	Fax 1-800-632-2109

Product Code	Size (cm)	Nomnal Thickness (mm)	Approximate water permeability* ml cm-2 min-1	Quantity per box
007937	2.5 x 10.2	1.50	3800 max	10
007939	5.1 x 10.2	1.50	3800 max	10
007826	10.2 x 10.2	1.50	3800 max	5
007827	15.2 x 15.2	1.50	3800 max	5

* Water Permeability data is approximate, data on file 1998.

† Bard and DeBakey are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.