



## Bard† Low Porosity PTFE Felt

### Ordering Information for Canada

WESTERN CANADA	Tel. 800 268 2862	Fax 1-800-632-2109
ONTARIO	Tel. 800 387 9473	Fax 1-800-632-2109
EASTERN CANADA	Tel. 800 387 7851	Fax 1-800-632-2109

<b>Product Code</b>	<b>Size (cm)</b>	<b>Nominal Thickness (mm)</b>	<b>Approximate water permeability* ml cm-2 min-1</b>	<b>Quantity per box</b>
007021	2.5 x 2.5	1.85	< 550 mean	5
007838	10.2 x 10.2	1.85	< 550 mean	5
007839	15.2 x 15.2	1.85	< 550 mean	5

\* Water Permeability data is approximate, data on file 1998.

† Bard is a registered trademark of C.R. Bard, Inc.