



## Ultraverse† Small Vessel PTA Catheters

### Ordering Information for Canada

WESTERN CANADA	Tel. 800 268 2862	Fax 1-800-632-2109
ONTARIO	Tel. 800 387 9473	Fax 1-800-632-2109
EASTERN CANADA	Tel. 800 387 7851	Fax 1-800-632-2109

75 cm	Shaft Length		Balloon Diameter (mm)	Size Length (cm)	RBP* (atm)	Nominal Pressure (atm)	Sheath Size
	100 cm	120 cm					
	UV-10022	UV-12022	2	2	15	8	4F
	UV-10024	UV-12024	2	4	15	8	4F
	UV-10026	UV-12026	2	6	15	8	4F
	UV-100210	UV-120210	2	10	10	8	5F
	UV-10032	UV-12032	3	2	15	8	5F
	UV-10034	UV-12034	3	4	15	8	4F
	UV-10036	UV-12036	3	6	15	8	4F
	UV-100310	UV-120310	3	10	10	8	5F
	UV-10042	UV-12042	4	2	12	8	4F
	UV-10044	UV-12044	4	4	12	8	4F
	UV-10046	UV-12046	4	6	12	8	4F
	UV-100410	UV-120410	4	10	10	8	5F
	UV-10052	UV-12052	5	2	12	8	5F
	UV-10054	UV-12054	5	4	10	8	5F
	UV-10056	UV-12056	5	6	10	8	5F
	UV-100510	UV-120510	5	10	10	8	5F
UV-7562		UV-12062	6	2	12	8	5F
UV-7564		UV-12064	6	4	12	8	5F
UV-7572		UV-12072	7	2	12	8	5F
UV-7574		UV-12074	7	4	12	8	5F

\* RBP (Rated Burst Pressure): the pressure at which BARD has 95% confidence that 99.9% of the balloons will not burst at or below upon single inflation.

\*\* Please contact your Bard Sales Representative for availability of sizes.

† Bard and Ultraverse are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.